

ROCCIA[®] ACIF ANTERIOR CERVICAL INTERBODY FUSION

INSTRUMENTATION GUIDE



MADE IN GERMANY

TABLE OF CONTENTS

Preface.....	3
Product information.....	4
INSTRUMENTATION	5
Position and approach; distraction	6
Discectomy and decompression	7
Selecting the trial implant	8
Filling the ROCCIA ACIF Cage	10
Inserting the ROCCIA ACIF Cage	11
Removing the instrument set.....	12
PRODUCT INFORMATION.....	14
Implants.....	PI 02 + 04
Trial implants	PI 03 + 05
Instruments	PI 06

NOTE: This Guide describes the use of ROCCIA® ACIF for anterior cervical instrumentation – it does not replace briefing by a surgeon experienced in surgical instrumentation of the spinal column.

We would be happy to assist you in finding a hospital that provides an opportunity to observe surgical procedures.

PREFACE

ROCCIA[®] ACIF – ANTERIOR CERVICAL INTERBODY FUSION

ROCCIA ACIF was developed for primary stabilization and restoration of the physiological lordosis in the cervical spine. The cage is designed for anterior approaches.

The chamber system in the cage improves intercorporal fusion as its generously proportioned design allows for the insertion of either bone or bone graft materials. At the same time, the cage has a broad supporting surface that largely prevents sinking when implanted correctly.

The broad portfolio with different heights, supporting surfaces and shapes allows individual selection of the implant based on the patient's anatomy.

Pins anchored in the cage allow good visualization of the correct position of the implant in the image converter and aid with the primary stability.

An intuitive instrument set with only a few surgical steps enables easy and efficient insertion of the cage.



NOTE: ROCCIA ACIF can be combined with additional stabilization depending on the stability and the sagittal profile.

NOTE: Please also follow the Instructions for Use provided with each product. All instrumentation guides and Instructions for Use can be found on our eLabeling portal (<https://elabeling.silony-medical.com/>).

ROCCIA® ACIF

PROVEN SYSTEM – COMBINED AND OPTIMISED

ROCCIA ACIF – The reliable cage for individual treatment




Preparing the fusion:

There is plenty of room in the generously proportioned ROCCIA ACIF chamber system for autologous or homologous bone and/or bone graft material for subsequent fusion.



Save time:

The delicate two-part ROCCIA ACIF Inserter is particularly easy to clean and integrates optimally into the lean instrument set. The black coating allows the inserter to be used under the microscope without producing reflections.




Rely on safety:

The ROCCIA ACIF Depth Stop guarantees you safe and targeted insertion of the cage and provides a good option to control the position of the cage.



For firm anchorage:

Symmetrical tapered pyramidal teeth on the surface of the PEEK implant along with radiographic markers ensure good primary anchoring of the cage to the base and cover plates.



Keep control:

Radiographic markers facilitate imaging of the ACIF implant during insertion as well as for the postoperative follow-up examination.



For individual treatment:

The implant portfolio offers a wide range of sizes and shapes (5°, wedge-shaped & anatomic), enabling you to provide individual treatment depending on the patient's anatomy.

ROCCIA[®] ACIF INSTRUMENTATION

The following section describes each of the necessary steps when using the ROCCIA ACIF for anterior cervical interbody fusion. Using ROCCIA ACIF, both monosegmental and multisegmental treatments can be carried out in accordance with this instrumentation guide.

Instrumentation of ACIF – Position and approach; distraction

The patient is placed in a supine position. The head should be stably positioned in a slightly reclined position. When positioning the patient, care should be taken to ensure that the target segment can be depicted well in the fluoroscopy both laterally and in A-P projection.

The approach is carried out with the customary procedure used in anterior cervical surgery. Standard retractors (e.g. Caspar retractor) support the direct and complete exposure of the target segment (Figs. 1 and 2).

In order to facilitate the implantation of the ROCCIA ACIF Cage, good exposure and distraction is recommended. To accomplish this, superior and inferior distraction pins should be set in parallel to the corresponding end plates of the target segment. The distractor can then be mounted.

Distractors of this type have supportive functions:

- Distraction of the target segment
- Stability of the target segment throughout the entire surgical procedure
- Parallel alignment of the vertebrae




Fig. 1 Exposure of the target segment using standard retractors




Fig. 2 Exposure of the target segment by spreading the retractor

Discectomy and decompression

The cervical intervertebral disc must be removed completely (Fig. 3). Therefore, box-shaped resection of the anterior longitudinal ligament is necessary first. Discectomy is performed as usual.

The cartilages on the end plates should be removed thoroughly until the end plates start to bleed slightly. When doing this, care should be taken to ensure that the end plates are not weakened in order to ensure sufficient supporting surface and stability for the implantation of the ROCCIA ACIF Cage.

The neural structures are decompressed using punches or high-speed milling cutters. In order to achieve good access to the target segment for the implantation, anterior osteophytes or other bony changes must also be removed if necessary.




Fig. 3 Complete removal of the cervical intervertebral disc

NOTE: Careful preparation of the disc space, especially extensive cleaning of the end plates, provides the basis for better vascularization and successful bony fusion. Damage of the bony base and cover plate can lead to sinking of the implant into the vertebral body.

Selecting the trial implant

RI-8010*
ROCCIA ACIF Inserter,
dismountable



RI-8050
ROCCIA ACIF Depth Stop, lat.
large



RI-T06161351**
ROCCIA ACIF Trial 6×16×13 mm,
anatomic



The ROCCIA ACIF Cage System offers a broad selection of different lengths, widths, angles and anatomic shapes, each with 5° lordosis. This extensive portfolio enables individual customization to different patient anatomies and intra-operative requirements.

With the ROCCIA ACIF Trial Implant, you can determine the implant to be used according to the individual anatomic situation. With the aid of the trial implant, you define the length, width and height and at the same time check which anatomic shape is suitable for the situation.

First, the insertion instrument must be assembled (Fig. 4). To do this, the shaft is inserted completely into the sleeve and screwed on tightly. The depth stop is then attached onto the instrument (Fig. 4). The depth stop can be removed and re-attached without removing the trial implant

* Representative for other inserters

** Representative for other trial implant sizes
see ROCCIA Trial Implants




Fig. 4 Assembly of the insertion instrument, consisting of the inserter, shaft and depth stop

Selecting the trial implant

RI-8010*
ROCCIA ACIF Inserter,
dismountable



RI-8050
ROCCIA ACIF Depth Stop, lat.
large



RI-T06161351**
ROCCIA ACIF Trial 6×16×13 mm,
anatomic



The desired trial implant is screwed onto the inserter (Fig. 5) and inserted into the intervertebral space under lateral X- ray control (Figs. 6 and 7). If anatomically shaped trials are used, the convex end plate must always face toward the superior direction. The label 'TOP' on the inserter must face the superior direction when inserting the trial (Fig. 6 zoom).

Silony Spine recommends using an implant that is as wide as possible to achieve a large contact surface area and to ensure support on the anterior and posterior cortical region of the end plates.

To determine the height, it is important to make sure that the trial implant is neither too tight nor too loose. Accordingly, try out a smaller or larger version until the trial implant is seated stably in the intersegmental space. When doing this, you should, if necessary, undo the distractor if you are using one, in order to get a tactile feel for the primary stability of the implant.

* Representative for other inserters

** Representative for other trial implant sizes
see ROCCIA Trial Implants

NOTE: The external dimension of the trial implant corresponds to the core dimension of the implant without the interlock or the X-ray marker.

NOTE: The trial implant also serves to simulate the insertion of the implant into its final position. If the trial implant cannot be brought into its desired end position, then renewed preparation of the intervertebral disc space may be necessary.




Fig. 5 Fixing the trial implant onto the inserter




Fig. 6 Introducing the trial implant




Fig. 7 X-ray control during insertion of the trial implant

Filling the ROCCIA ACIF

RI-2060
ROCCIA ACIF Loading



RI-2061
ROCCIA ACIF Bone Graft



RI-8010*
ROCCIA ACIF
Inserter.




Fig. 8 Fixing the cage onto the

The selected ROCCIA ACIF Cage is screwed onto the previously assembled ROCCIA Inserter (Fig. 8). The inserter must not be screwed too tightly onto the cage. When doing this, the depth stop can be left on the instrument.

Filling the cage with autologous bone material (or homologous bone or bone graft material) is an important precondition for reliable fusion. A loading block and a pusher are provided for this (Fig. 9). The color rings on the pusher correspond to the color of the previously used trial im-plant. The pusher surface is adjusted to the corresponding cage footprint. The alignment of the pusher corresponding to the graduation markings on the pusher and on the loading block must be noted for this purpose (Fig. 9 zoom). The graduation marking indicates the posterior side of the cage.

* Representative for other inserters
see ROCCIA instruments




Fig. 9 Filling the cage with bone material in the loading block with pusher

Inserting the ROCCIA ACIF Cage

RI-8010*
ROCCIA ACIF Inserter,
dismountable



RI-8050
ROCCIA ACIF Depth Stop, lat.
large



RI-8060
ROCCIA Driving Mallet, small



Check the superior and inferior alignment of the cage. If anatomically shaped cages are used, the convex end plate must always face the superior direction. The label 'TOP' on the inserter must face the superior direction when inserting the cage (Fig. 10 zoom).

The cage is inserted into the disc space under X-ray control (Fig. 10).

The depth stop on the inserter prevents the implant being inserted too deeply.

A small driving mallet is available to make insertion of the cage easier.

* Representative for other inserters
see ROCCIA instruments




Fig. 10 Inserting the cage

NOTE: The trial implant has a smooth surface. Because of the serration used for anchorage, the cage has a rough surface. This means that greater force may be required to insert the cage. If the surgical intervention is carried out with the aid of a distractor, then it is advisable to slightly increase the distraction in order to reduce the forces required for inserting the cage. After complete implantation, this distraction should immediately be loosened again.

NOTE: Correctly selecting the cage size has a decisive impact on the success of the instrumentation and fusion.

Removing the instrument set

The final position of the implant (Fig. 11) should be checked using the fluoroscopy (lateral and anterior-posterior view). When doing this, X-ray markers in the implant show the position of the cage.

After the final position has been confirmed in the fluoroscopy, the inserter is unscrewed completely from the implant and removed (Fig. 12 and 13). No force needs to be exerted to withdraw the inserter from the implant. Should resistance nevertheless be felt, it must be checked whether the inserter has been unscrewed completely. If necessary, this process must be repeated. Make sure that the final position of the implant is not altered when the inserter is removed.




Fig. 11 Final position of the cage




Fig. 12 X-ray control, lateral




Fig. 13 X-ray control, anterior-posterior


NOTE: Depending on the patient's individual pathology, additional stabilization of the ROCCIA ACIF Cage with an anterior plate or a posterior screw-rod system may be required.

ROCCIA® ACIF IMPLANTS PRODUCT INFORMATION

ROCCIA ACIF Anatomic implants by article number.....	PI 02
ROCCIA ACIF Anatomic trial implants by article number.....	PI 03
ROCCIA ACIF Wedge-shaped Implants by article number.....	PI 04
ROCCIA ACIF Wedge-shaped trial implants by article number.....	PI 05
ROCCIA ACIF Instruments by article number.....	PI 06

ROCCIA® ACIF Implants

Article number explanation for the cages, as an example.



System:
ROCCIA

Implant type:
ACIF

Configuration:
anatomic

Material:
PEEK
Ti6Al4V ELI (x-ray marker)



Article number	Description	Illustration
S-RCP-04141351-S	ROCCIA ACIF Cage, 4×14×13 mm, anatomic	
S-RCP-05141351-S	ROCCIA ACIF Cage, 5×14×13 mm, anatomic	
S-RCP-06141351-S	ROCCIA ACIF Cage, 6×14×13 mm, anatomic	
S-RCP-07141351-S	ROCCIA ACIF Cage, 7×14×13 mm, anatomic	
S-RCP-08141351-S	ROCCIA ACIF Cage, 8×14×13 mm, anatomic	
S-RCP-04161351-S	ROCCIA ACIF Cage, 4×16×13 mm, anatomic	
S-RCP-05161351-S	ROCCIA ACIF Cage, 5×16×13 mm, anatomic	
S-RCP-06161351-S	ROCCIA ACIF Cage, 6×16×13 mm, anatomic	
S-RCP-07161351-S	ROCCIA ACIF Cage, 7×16×13 mm, anatomic	
S-RCP-08161351-S	ROCCIA ACIF Cage, 8×16×13 mm, anatomic	
S-RCP-04161551-S	ROCCIA ACIF Cage, 4×16×15 mm, anatomic	
S-RCP-05161551-S	ROCCIA ACIF Cage, 5×16×15 mm, anatomic	
S-RCP-06161551-S	ROCCIA ACIF Cage, 6×16×15 mm, anatomic	
S-RCP-07161551-S	ROCCIA ACIF Cage, 7×16×15 mm, anatomic	
S-RCP-08161551-S	ROCCIA ACIF Cage, 8×16×15 mm, anatomic	
S-RCP-04181351-S	ROCCIA ACIF Cage, 4×18×13 mm, anatomic	
S-RCP-05181351-S	ROCCIA ACIF Cage, 5×18×13 mm, anatomic	
S-RCP-06181351-S	ROCCIA ACIF Cage, 6×18×13 mm, anatomic	
S-RCP-07181351-S	ROCCIA ACIF Cage, 7×18×13 mm, anatomic	
S-RCP-08181351-S	ROCCIA ACIF Cage, 8×18×13 mm, anatomic	
S-RCP-04181551-S	ROCCIA ACIF Cage, 4×18×15 mm, anatomic	
S-RCP-05181551-S	ROCCIA ACIF Cage, 5×18×15 mm, anatomic	
S-RCP-06181551-S	ROCCIA ACIF Cage, 6×18×15 mm, anatomic	
S-RCP-07181551-S	ROCCIA ACIF Cage, 7×18×15 mm, anatomic	
S-RCP-08181551-S	ROCCIA ACIF Cage, 8×18×15 mm, anatomic	

ROCCIA® ACIF Trial Implants

Illustration	Description	Article number
	ROCCIA ACIF Trial 4×14×13 mm, anatomic	RI-T04141351
	ROCCIA ACIF Trial 5×14×13 mm, anatomic	RI-T05141351
	ROCCIA ACIF Trial 6×14×13 mm, anatomic	RI-T06141351
	ROCCIA ACIF Trial 7×14×13 mm, anatomic	RI-T07141351
	ROCCIA ACIF Trial 8×14×13 mm, anatomic	RI-T08141351
	ROCCIA ACIF Trial 4×16×13 mm, anatomic	RI-T04161351
	ROCCIA ACIF Trial 5×16×13 mm, anatomic	RI-T05161351
	ROCCIA ACIF Trial 6×16×13 mm, anatomic	RI-T06161351
	ROCCIA ACIF Trial 7×16×13 mm, anatomic	RI-T07161351
	ROCCIA ACIF Trial 8×16×13 mm, anatomic	RI-T08161351
	ROCCIA ACIF Trial 4×16×15 mm, anatomic	RI-T04161551
	ROCCIA ACIF Trial 5×16×15 mm, anatomic	RI-T05161551
	ROCCIA ACIF Trial 6×16×15 mm, anatomic	RI-T06161551
	ROCCIA ACIF Trial 7×16×15 mm, anatomic	RI-T07161551
	ROCCIA ACIF Trial 8×16×15 mm, anatomic	RI-T08161551
	ROCCIA ACIF Trial 4×18×13 mm, anatomic	RI-T04181351
	ROCCIA ACIF Trial 5×18×13 mm, anatomic	RI-T05181351
	ROCCIA ACIF Trial 6×18×13 mm, anatomic	RI-T06181351
	ROCCIA ACIF Trial 7×18×13 mm, anatomic	RI-T07181351
	ROCCIA ACIF Trial 8×18×13 mm, anatomic	RI-T08181351
	ROCCIA ACIF Trial 4×18×15 mm, anatomic	RI-T04181551
	ROCCIA ACIF Trial 5×18×15 mm, anatomic	RI-T05181551
	ROCCIA ACIF Trial 6×18×15 mm, anatomic	RI-T06181551
	ROCCIA ACIF Trial 7×18×15 mm, anatomic	RI-T07181551
	ROCCIA ACIF Trial 8×18×15 mm, anatomic	RI-T08181551

System:
ROCCIA

Instrument
type: Trial
implant



Configuration:
anatomic

Material:
Ti6Al4V ELI

 D30003

ROCCIA® ACIF Implants

Article number explanation for the Cage as an example.




System:
ROCCIA

Implant type:
ACIF

Configuration:
wedge-shaped

Material:
PEEK
Ti6Al4V ELI (x-ray marker)









Article number	Description	Illustration
S-RCP-04141305-S	ROCCIA ACIF Cage, 4×14×13 mm, wedge-shaped	
S-RCP-05141305-S	ROCCIA ACIF Cage, 5×14×13 mm, wedge-shaped	
S-RCP-06141305-S	ROCCIA ACIF Cage, 6×14×13 mm, wedge-shaped	
S-RCP-07141305-S	ROCCIA ACIF Cage, 7×14×13 mm, wedge-shaped	
S-RCP-08141305-S	ROCCIA ACIF Cage, 8×14×13 mm, wedge-shaped	
S-RCP-04161305-S	ROCCIA ACIF Cage, 4×16×13 mm, wedge-shaped	
S-RCP-05161305-S	ROCCIA ACIF Cage, 5×16×13 mm, wedge-shaped	
S-RCP-06161305-S	ROCCIA ACIF Cage, 6×16×13 mm, wedge-shaped	
S-RCP-07161305-S	ROCCIA ACIF Cage, 7×16×13 mm, wedge-shaped	
S-RCP-08161305-S	ROCCIA ACIF Cage, 8×16×13 mm, wedge-shaped	
S-RCP-04161505-S	ROCCIA ACIF Cage, 4×16×15 mm, wedge-shaped	
S-RCP-05161505-S	ROCCIA ACIF Cage, 5×16×15 mm, wedge-shaped	
S-RCP-06161505-S	ROCCIA ACIF Cage, 6×16×15 mm, wedge-shaped	
S-RCP-07161505-S	ROCCIA ACIF Cage, 7×16×15 mm, wedge-shaped	
S-RCP-08161505-S	ROCCIA ACIF Cage, 8×16×15 mm, wedge-shaped	
S-RCP-04181305-S	ROCCIA ACIF Cage, 4×18×13 mm, wedge-shaped	
S-RCP-05181305-S	ROCCIA ACIF Cage, 5×18×13 mm, wedge-shaped	
S-RCP-06181305-S	ROCCIA ACIF Cage, 6×18×13 mm, wedge-shaped	
S-RCP-07181305-S	ROCCIA ACIF Cage, 7×18×13 mm, wedge-shaped	
S-RCP-08181305-S	ROCCIA ACIF Cage, 8×18×13 mm, wedge-shaped	
S-RCP-04181505-S	ROCCIA ACIF Cage, 4×18×15 mm, wedge-shaped	
S-RCP-05181505-S	ROCCIA ACIF Cage, 5×18×15 mm, wedge-shaped	
S-RCP-06181505-S	ROCCIA ACIF Cage, 6×18×15 mm, wedge-shaped	
S-RCP-07181505-S	ROCCIA ACIF Cage, 7×18×15 mm, wedge-shaped	
S-RCP-08181505-S	ROCCIA ACIF Cage, 8×18×15 mm, wedge-shaped	

ROCCIA® ACIF Trial Implants

Illustration	Description	Article number	System: ROCCIA
	ROCCIA ACIF Trial, 4×14×13 mm, wedge-shaped	RI-T04141305	Instrument type: Trial implant
	ROCCIA ACIF Trial, 5×14×13 mm, wedge-shaped	RI-T05141305	
	ROCCIA ACIF Trial, 6×14×13 mm, wedge-shaped	RI-T06141305	
	ROCCIA ACIF Trial, 7×14×13 mm, wedge-shaped	RI-T07141305	
	ROCCIA ACIF Trial, 8×14×13 mm, wedge-shaped	RI-T08141305	
	ROCCIA ACIF Trial, 4×16×13 mm, wedge-shaped	RI-T04161305	Configuration: wedge- shaped
	ROCCIA ACIF Trial, 5×16×13 mm, wedge-shaped	RI-T05161305	
	ROCCIA ACIF Trial, 6×16×13 mm, wedge-shaped	RI-T06161305	
	ROCCIA ACIF Trial, 7×16×13 mm, wedge-shaped	RI-T07161305	
	ROCCIA ACIF Trial, 8×16×13 mm, wedge-shaped	RI-T08161305	
	ROCCIA ACIF Trial, 4×16×15 mm, wedge-shaped	RI-T04161505	Material: Ti6Al4V ELI
	ROCCIA ACIF Trial, 5×16×15 mm, wedge-shaped	RI-T05161505	
	ROCCIA ACIF Trial, 6×16×15 mm, wedge-shaped	RI-T06161505	
	ROCCIA ACIF Trial, 7×16×15 mm, wedge-shaped	RI-T07161505	
	ROCCIA ACIF Trial, 8×16×15 mm, wedge-shaped	RI-T08161505	
	ROCCIA ACIF Trial, 4×18×13 mm, wedge-shaped	RI-T04181305	⚠ D30003
	ROCCIA ACIF Trial, 5×18×13 mm, wedge-shaped	RI-T05181305	
	ROCCIA ACIF Trial, 6×18×13 mm, wedge-shaped	RI-T06181305	
	ROCCIA ACIF Trial, 7×18×13 mm, wedge-shaped	RI-T07181305	
	ROCCIA ACIF Trial, 8×18×13 mm, wedge-shaped	RI-T08181305	
	ROCCIA ACIF Trial, 4×18×15 mm, wedge-shaped	RI-T04181505	
	ROCCIA ACIF Trial, 5×18×15 mm, wedge-shaped	RI-T05181505	
	ROCCIA ACIF Trial, 6×18×15 mm, wedge-shaped	RI-T06181505	
	ROCCIA ACIF Trial, 7×18×15 mm, wedge-shaped	RI-T07181505	
	ROCCIA ACIF Trial, 8×18×15 mm, wedge-shaped	RI-T08181505	

ROCCIA® ACIF Instruments

Article number	Description	Illustration
RI-2060*	ROCCIA ACIF Loading Block	
RI-2061*	ROCCIA ACIF Bone Graft Pusher	
RI-8010	ROCCIA ACIF Inserter, dismountable	
RI-8050	ROCCIA ACIF Depth Stop, lat. large	
RI-8060*	ROCCIA Driving Mallet, small	
RI-8110	ROCCIA ACIF Inserter, long, dismountable	



 Silony Medical GmbH
Leinfelder Straße 60
70771 Leinfelden-Echterdingen
Germany
Tel +49 711 78 25 25 0
Fax +49 711 78 25 25 11

 elabeling.silony-medical.com/contact

www.silonyspine.com

D30038.d.EN 12.06.2025

CE 0483

