(en)
INSTRUCTION FOR USE (IFU)
Modular Retractor System (MRS)

  33.430.205
0476
MEDICALEX
94110 Arcueil
France
+33 (0)1 46 11 16 20
www.medicalcx.com/info

  10-306-16 (Imech 5067M288)
Reison Medical AB
Erikbergsvägen 32A
734 92 Hallstahammar
Sweden
<http://reison.se/se>

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Modular Retractor System (MRS)

INTRODUCTION:

As Bradshaw Medical Inc. (dba Intech), we want you to have the best possible experience with our instruments and to be able to use them with the highest standards in the medical industry. We have designed and produced this kit with the highest standards, and to make it intuitive to use. Therefore, please read this guide for your use. This document was created to provide the user with detailed instructions on how to assemble the system together, as well as how to maintain the device and determine when maintenance is needed.

Bradshaw Medical Inc. (dba Intech) manufactures reusable devices for Spine and Orthopedic surgical activities. The Modular Retractor System (MRS) is a set of reusable surgical instruments designed for transient use in orthopedic surgical procedures. To ensure the device has this functionality, the MRS consists of multiple components listed below in Devices in scope section that will be assembled and used by a trained orthopedic surgeon or neurosurgeon to create a posterior access to the spine.

INTENDED USE:

The Modular Retractor System (MRS) is designed to perform surgery on patients with spine pathology. The MRS is a medical device system designed to distract and hold vertebrae, muscle and/ or soft tissue to enable passage of instruments that will be used during spinal surgery. It is designed for the needs of various indications and/or surgical techniques.

The following devices and accessories intended use are listed below:

Retractor: This is a medical device designed to separate and hold soft tissue/muscle and vertebrae dissection in order to enable passage of instrumentation during orthopedic and spinal MIS (Minimal Invasive Surgery) surgery. The Retractor is the Rack, Bridge, Distraction compression case, Fixed case, Mobile case, Articulated arm, Right rotative arm, and Left rotative arm.

Dilator: A surgical medical device used to induce dilation which expands an opening from a surgical incision.

Threaded Shim: A two-piece surgical retractor fixation device/shim for use in surgical procedures, particularly spinal surgeries. The surgical retractor fixation device/shim is designed to lock a retractor to an anatomical structure, such as the track of the vertebral body.

Tray: A medical device family with multi levels that holds all the surgical tools and equipment that is expected to be required to complete a surgical procedure. The surgical instruments are cleaned, placed in the tray for sterilization and to store and secure the device during shipping. The tray medical device family is Tray 2, Lid 1, Lid 2, Base 1, and Base 2.

Driver, 3.5mm hex: A screwdrivers is a surgical medical device specifically used to screw/unscREW the threaded shaft of the Threaded Shim onto the vertebrae.

Wrench, 5mm hex: The wrench is a surgical medical device specifically used to facilitate the adjustment and the motion of the retractor.

Concentrical Blades 50 mm to 130 mm (9 lengths): A surgical blade in a concentrical shape used to separate and hold soft tissue/muscle to enable passage of instruments that will be used during orthopedic/spinal Minimal Invasive Surgery (MIS).

Plate Blades w/ teeth 50 mm to 130mm (9 lengths): A flat surgical blade with teeth used to separate and hold soft tissue/muscle to enable passage of instruments that will be used during orthopedic/spinal Minimal Invasive Surgery (MIS).

Plate Blades without teeth 50 mm to 130mm (9 lengths): A flat surgical blade used to separate and hold soft tissue/muscle to enable passage of instruments that will be used during orthopedic/spinal Minimal Invasive Surgery (MIS).

PATIENT GROUP:
Due to the nature of this device, there are no specific demographic that is assigned to this device.

INDICATIONS and INTENDED PURPOSE:

The MRS kit does not have an indicative characteristic. This device is used to facilitate the required surgical intervention. Bradshaw Medical Inc. (dba Intech) is not aware of any contraindications with regards to the MRS kit. In cases where the MRS is used in combination with implants or instruments, please refer to the respective instructions for use of the related devices and additional surgical steps.

CONTRAINDICATIONS:

The MRS does not have an indicative characteristic. This device is used to facilitate the required surgical intervention. Bradshaw Medical Inc. (dba Intech) is not aware of any contraindications with regards to MRS but continues to monitor the use and application of these devices. Use as intended only. Don't use the MRS kit for surgeries other than posterior spinal surgery. Use only the wrenches available in the MRS kit for the motion of the MRS to prevent overload.

INTENDED USER:


Surgical procedures should be performed only by experienced spinal surgeons having adequate training and familiarities with surgical techniques. Medical staff must be trained for proper assembly of the system and all instruction regarding safety and usage must be read carefully before use.

DEVICE DESCRIPTION:

The surgical set is comprised of components required to make the needed retraction device, using the usage instructions provided below for Posterior Spinal Surgery. Devices are supplied **NON-STERILE** and must be inspected, cleaned, and sterilized before each use.

DEVICES IN SCOPE:

Illustration of the fully loaded set:


Medical Device Part Number	Description	Quantity
MRS-100	Modular Retractor System	1
MRS-101	Rack	2
MRS-102	Bridge	1
MRS-103	Distraction compression case	2
MRS-104-01	Fixed case	2
MRS-104-02	Mobile case	2
MRS-105	Articulated arm	2
MRS-106-R	Right rotative arm	1
MRS-106-L	Left rotative arm	1
MRS-107-1	Dilator 1	1
MRS-107-2	Dilator 2	1
MRS-107-3	Dilator 3	1
MRS-111	Threaded shim	2
MRS-114	Wrench, 5mm hex	1
MRS-116-01-01	Lid 1	1
MRS-116-01-02	Base 1	1
MRS-116-02-01	Lid 2	1
MRS-116-02-02	Tray 2	1
MRS-116-02-03	Base 2	1
MRS-118	Driver, 3.5mm hex	1

Medical Device Part Number	Description	Quantity
MRS-108-050	Concentrical blade - length 50mm	2
MRS-108-060	Concentrical blade - length 60mm	2
MRS-108-070	Concentrical blade - length 70mm	2
MRS-108-080	Concentrical blade - length 80mm	2
MRS-108-090	Concentrical blade - length 90mm	2
MRS-108-100	Concentrical blade - length 100mm	2
MRS-108-110	Concentrical blade - length 110mm	2
MRS-108-120	Concentrical blade - length 120mm	2
MRS-108-130	Concentrical blade - length 130mm	2
MRS-109-050	Plate blade with teeth - length 50 mm	2
MRS-109-060	Plate blade with teeth - length 60 mm	2
MRS-109-070	Plate blade with teeth - length 70 mm	2
MRS-109-080	Plate blade with teeth - length 80 mm	2
MRS-109-090	Plate blade with teeth - length 90 mm	2
MRS-109-100	Plate blade with teeth - length 100 mm	2
MRS-109-110	Plate blade with teeth - length 110 mm	2
MRS-109-120	Plate blade with teeth - length 120 mm	2
MRS-109-130	Plate blade with teeth - length 130 mm	2
MRS-110-050	Plate blade without teeth - length 50 mm	2
MRS-110-060	Plate blade without teeth - length 60 mm	2
MRS-110-070	Plate blade without teeth - length 70 mm	2
MRS-110-080	Plate blade without teeth - length 80 mm	2
MRS-110-090	Plate blade without teeth - length 90 mm	2
MRS-110-100	Plate blade without teeth - length 100 mm	2
MRS-110-110	Plate blade without teeth - length 110 mm	2
MRS-110-120	Plate blade without teeth - length 120 mm	2
MRS-110-130	Plate blade without teeth - length 130 mm	2

Medical Device Accessory Part Number	Description
S061TM293	Retractor illuminator lightguide-bifurcated
S061TM294	Adapter, bifurcated ACMI-UNIVERSAL F.O. cable
MRS-113	INTECH 40mm surgical illuminator
S061TM234	Adapter, bifurcated ACMI F.O. cable
S061TM278	INTECH STORZ-UNIV adaptor
S061TM268	INTECH WOLF-UNIV adaptor
S061TM269	INTECH OLYMPUS-UNIV adaptor
S061TM270	INTECH OLYMPUS-UNIV adaptor
S061TM288	Radial Clamp (Table Clamp)
33-430-205	Broche Kirchner 1. pte Trocard 30cm 2mm
MRS-115	Arm (Articulated Arm for Radial Clamp)

Usage Instructions:


Connection :


- Align the arm with the rack's quick connect and insert it into the Rack connection hole until you hear a click

- Press on the Release Button while pulling on the connected part (rotative arm for illustration) for removal

Connection: to be connected to the rack MRS-101 or compression distraction case MRS-103

Selector Mode

Turn the selector to switch mode from Compression /lock / distraction as shown below:

Connection of the rack

Choose the appropriate position on the mobile case (compression or distraction) to insert the mobile case onto the rack


- Select the lock position to immobilize the mobile case

Several mobile case can be connected onto the rack

Connection:



- Align the arm with the rack's quick connect and insert it into the Rack connection hole until a click is heard

Disconnection:


- Press on the Release Button while pulling on the connected part (relative arm for illustration) for removal

Actuation of the Compression Distraction Case:

Insert Wrench into the Hex footprint to move the case in compression or distraction

Connection of the mobile case to the fixed case:


- Align rack with the housing to insert it into the fixed case (refer to the laser marked shape for the orientation)



- Push the release button to disassemble while pulling on the rack . . .

Sliding the mobile case along the rack:


Connection of the fixed case to the rack:

- Choose the position of the case onto the rack
- Simply insert the fixed case onto the rack until you hear a click
- Push on the Release Button to move the case on the rack or for disassembly


Connection of a blade to the mobile case:

- Insert blade into arm until you hear a click

- Push on the Release Button and pull the blade out for disassembly

Tilting the blade:

Articulated Arm MRS-105


Connection: to be connected to the rack or compression distraction case

Tips: to be connected to a screw


Connection of the rack or compression distraction case:

- Refer to MRS-101 "connecting / disconnecting a part"

Movement of the arm:

Right/Left Rotative Arm MRS-106-R/L


Locking button: to lock and unlock the articulation of the arm
Release button: to unlock the connected blade

Connection: to be connected to the rack or compression distraction case
Double hex imprint: to tilt the tip

Connecting to the rack or compression distraction case:


- Refer to MRS-101 "connecting / disconnecting a part"

Connecting a blade to the rotative arm:


- Present the Arm MRS-115 in front of the hole
- Push the part until you hear the click
- Push the release button
- Pull on the blade to disassemble

Tilting of the arm:

- Put the wrench in the Double Hex imprint
- Turn the wrench to tilt the arm

Rotation of the arm:


- Pull the Lock Button
- Move the arm
- To immobilize the arm, put the arm straight then push the lock button

Dilator MRS-107-01/02/03


Depth indicator: to indicate the depth of the access

Assembling of the dilator

- Assembling the dilator from the smallest to the biggest one.
- A K-Wire can be used as a guide.

Concentrical Blade MRS-108

Hole for light: to attach the light
Dove-tail: to attach the illuminator lightguide-bifurcated or the shim

Blade attachment: to attach the blade to the retractor

Use with the dilators

- Use 2 identical concentric blades

- Slide around the dilators

Connection of the blade:

- Refer to "Fixed case MRS-104-01 + Mobile case MRS-104-02" or "Right/Left relative arm MRS-106-R/L"

Connection of the light cable S06ITW293:

- 
- Present the light tip in front of the hole
 - Slide the tip into the hole

Connection of the illuminator lightguide-bifurcated MRS-113:

- 
- Present the light mat on top of the blade
 - Slide down the tip into the dovetail

Connection of the Threaded Shim MRS-111:



- 
- Present the threaded shim on top of the blade
 - Slide down into the dovetail

Plate blade with teeth MRS-109


Connecting the light cable S06ITM293:

- Refer to MRS-108 “Connecting of the light cable S06ITM293”

Connection of the light cable S06ITM293:


- Refer to MRS-108 “Connecting of the light cable S06ITM293”

Screwing into a bone

- Use the screwdriver MRS-118 and the blade MRS-108
- Turn the handle to screw the threaded part into the vertebrae
- Refer to the blade MRS-108 for the shim connection


Connecting to the parts and motion:

Refer to:


- Compression Distraction case MRS-103 motion of the Compression Distraction case
- Fixed case MRS-104-01 + Mobile case MRS-104-02
- Right/Left rotative arm MRS-106-R/L

Driver: 3.5mm hex MRS-118


Connecting to the shim

- Refer to threaded shim MRS-111

Connecting the rack:

Present the arm MRS-115 in front of the hole

Introduce the pin in the hole

Ensure to first unscrew the locking knob (counterclockwise)

You can choose between three positions

Screw the locking knob for locking (clockwise)

Unscrew the locking knob for unlocking (counterclockwise)

Setting of the arm:

- Maintain the arm before unscrewing the central knob
- Unscrew the central knob for unlocking (counterclockwise)
- Set the arm in the position as needed
- Screw the locking knob for locking (clockwise)

Cleaning and Reprocessing Information:

<p>WARNINGS:</p> <ul style="list-style-type: none"> • Avoid contact with: Chemicals such as Hydrochloric Acid, Sulfuric Acid, Hydrofluoric Acid, Alkaline chemicals, and high concentrations of chlorides such as Chlorine Bleach. • No specific enzymatic detergent is endorsed but the final solution should have a neutral pH between 7.0 and 8.5 according to ASTM F1744-96, the standard guide for care and handling of stainless steel. • NEWLY PURCHASED DEVICES MUST UNDERGO DECONTAMINATION PRIOR TO STERILIZATION. • When using any detergent or equipment, always follow the manufacturers' instructions / recommendations for that product. 	<ul style="list-style-type: none"> • Avoid contact with: Chemicals such as Hydrochloric Acid, Sulfuric Acid, Hydrofluoric Acid, Alkaline chemicals, and high concentrations of chlorides such as Chlorine Bleach. • No specific enzymatic detergent is endorsed but the final solution should have a neutral pH between 7.0 and 8.5 according to ASTM F1744-96, the standard guide for care and handling of stainless steel. • NEWLY PURCHASED DEVICES MUST UNDERGO DECONTAMINATION PRIOR TO STERILIZATION. • When using any detergent or equipment, always follow the manufacturers' instructions / recommendations for that product.
<p>Limitations on Reprocessing: Repeated reprocessing has been included in the validation for the lifetime of the performance of this device.</p> <ul style="list-style-type: none"> • Clean after each surgery. • All devices should be removed from the tray. • Submerge the device(s) in an enzymatic detergent that is safe for medical devices. • Prepare the detergent according to manufacturer's recommendations. • Soak the device(s) for 5 minutes in the detergent. • Scrub the submerged device(s) with a soft bristled brush to remove any remaining soil or debris, paying special attention to small grooves, crevices, and cannulation. • Agitate the device(s) in the solution while scrubbing. Use a small cytology brush to clean internal channels and actuate any movable parts to loosen any trapped soil. Rinse thoroughly with warm (38-49°C / 100.4F-120.2F) demineralized water. • Agitate the device by hand for at least 3 minutes in a bath of warm (38-49°C / 100.4F-120.2F) demineralized water. • Check this device for any visible soil and repeat the cleaning process if soil is found. (A longer presoak may be necessary). • Ultrasonically clean the device(s) for 10 minutes in an enzymatic detergent (or acceptable alternative). • Prepare the detergent according to the manufacturer's recommendations. • Rinse the device(s) with clean demineralized water for at least 1 minute. If the water leaving the device appears to be soiled, repeat the cleaning procedure. • Dry the exterior of the device(s) with a clean, lint-free cloth. • DO NOT use any abrasive detergents, or cleaning pads. All brushes used should be designed for use on medical devices. 	<p>Limitations on Reprocessing: Repeated reprocessing has been included in the validation for the lifetime of the performance of this device.</p> <ul style="list-style-type: none"> • Clean after each surgery. • All devices should be removed from the tray. • Submerge the device(s) in an enzymatic detergent that is safe for medical devices. • Prepare the detergent according to manufacturer's recommendations. • Soak the device(s) for 5 minutes in the detergent. • Scrub the submerged device(s) with a soft bristled brush to remove any remaining soil or debris, paying special attention to small grooves, crevices, and cannulation. • Agitate the device(s) in the solution while scrubbing. Use a small cytology brush to clean internal channels and actuate any movable parts to loosen any trapped soil. Rinse thoroughly with warm (38-49°C / 100.4F-120.2F) demineralized water. • Agitate the device by hand for at least 3 minutes in a bath of warm (38-49°C / 100.4F-120.2F) demineralized water. • Check this device for any visible soil and repeat the cleaning process if soil is found. (A longer presoak may be necessary). • Ultrasonically clean the device(s) for 10 minutes in an enzymatic detergent (or acceptable alternative). • Prepare the detergent according to the manufacturer's recommendations. • Rinse the device(s) with clean demineralized water for at least 1 minute. If the water leaving the device appears to be soiled, repeat the cleaning procedure. • Dry the exterior of the device(s) with a clean, lint-free cloth. • DO NOT use any abrasive detergents, or cleaning pads. All brushes used should be designed for use on medical devices.
<p>Cleaning (Manual):</p> <ul style="list-style-type: none"> • Ultrasonically clean the device(s) for 10 minutes in an enzymatic detergent (or acceptable alternative). • Prepare the detergent according to the manufacturer's recommendations. • Rinse the device(s) with clean demineralized water for at least 1 minute. If the water leaving the device appears to be soiled, repeat the cleaning procedure. • Dry the exterior of the device(s) with a clean, lint-free cloth. • DO NOT use any abrasive detergents, or cleaning pads. All brushes used should be designed for use on medical devices. 	<p>Cleaning (Automated):</p> <ul style="list-style-type: none"> • Clean after each surgery. • All devices should be removed from the tray. • Prepare the detergent according to the manufacturer's recommendation. • Ultrasonically clean the device(s) for 10 minutes in a neutral pH detergent (or acceptable alternative). • Run the automatic wash cycle with a temperature of 40°C and ultrasonic waves of 35 kHz for 10 minutes. These times and temperatures will vary depending on the equipment and

Important Precautions and Safety Instructions:

1. Read the instructions for Use in their entirety before operating the Modular Retractor System. Failure to read the manual could result in harmful effects to the user, patient, and/or the MIS device.
2. Inspect devices before every use. Do not use if the device or parts of the device are damaged. Only use if the Modular Retractor System is in proper functional condition.
3. Discontinue use immediately if device is damaged or in nonfunctional condition.
4. Do not use the Modular Retractor System for any purpose other than the distraction or compression of the vertebrae, muscle, or soft tissue, during orthopedic/spine surgeries.
5. Do not use the device if any of the following event(s) occurs:
 - a. Do not use on any item
 - b. Failure at cleaning or sterilization
 - c. Excessive force required for any actuation of the device
 - d. Unable to loosen or tighten fasteners,
 - e. Excessive wear on any item
 - f. Blades do not properly fit into articulated arm.
 - g. Slippage of blade holder on the rack under normal conditions.
6. Precaution use to prevent damage:
 - a. Do not impact, even gently, any of the provided devices
 - b. Do not expose the devices and accessories to temperatures less than -10 °F (-23.3°C) or greater than 235 °F (162.8°C).
7. If there is any serious incident that occurs in relation to the devices, contact and report it to Bradshaw Medical (dba Intech Medical) and contact the competent authority and or the authorized representative (Emergo) of the market where the patient is located. **(See Return and Service Order Instructions)**


	<p>detergent use. Following cleaning, thoroughly rinse the device with demineralized water for at least 1 minute to remove any debris or detergents.</p> <ul style="list-style-type: none"> Repeat this cleaning procedure if the device(s) appear to be soiled after initial cleaning. Rinse the device(s) with clean tap water for at least 1 minute. If the water leaving the device appears to be soiled, repeat the cleaning procedure. Sterilization before each use 	<p>Gravity</p>
Recommended Sterilization	<p>Pre-Double wrapped</p> <p>Vacuum wrapped</p> <p>Cycle: 135°C (275°F)</p> <p>Temperature: 135°C (275°F)</p> <p>Exposure Time: 5 minutes</p> <p>Dry Time: 50 minutes</p> <p>Dry Time: 50 minutes</p> <p>*NOTE* After a steam sterilization cycle the device <u>must</u> be allowed to cool for a minimum of 40 minutes in (stagnant) open air prior to being used.</p>	<p>Double Wrapped</p> <p>135°C (275°F)</p> <p>15 minutes</p> <p>50 minutes</p>
Inspection:	<ul style="list-style-type: none"> Inspection after each cleaning and use. Ensure that there are no spots or stains on the device after cleaning/disinfection. If such spots exist, repeat the cleaning process to remove any protein residues or other stains. 	
Packaging:	<ul style="list-style-type: none"> Ensure the device has dried and is free of spots prior to packaging. Ensure that the packaging material has been held at room temperature for a minimum of 2 hours before use. Package the device using standard packaging materials and packaging practices accepted by the institution. 	
Storage:	<ul style="list-style-type: none"> The device should be stored in an area that protects them from dusts, moisture, insects, vermin, and extremes of temperature and humidity. Follow instruction's storage policies. 	

Return and Service Order Instructions:

Due to cases of misuse or device failure, the Modular Retractor System may require service or a return. If any wear, damage, or previously mentioned hazards are present, see below for information.

If the instrument is damaged for any reason, the user must return the device to the manufacturer for disposing of the device.

Symbols Glossary:

Reference Number of Symbol	Title of Symbol	Description of Symbol
	Manufacturer	Indicates the medical device manufacturer and is accompanied by the name and the address of the manufacturer
	Authorized representative in the European Community	Indicates the authorized representative in the European Community.
	Date of manufacture	Indicates the date when the medical device was manufactured.
	Batch code	Indicates the manufacturer's batch codes so that the batch or lot can be identified.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	Medical Device	Indicates that the product is the medical device
	Unique Device Identifier	Indicates a carrier that contains Unique Identifier Information
	Distributor	Indicates the entity distributing the medical device into the locale

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