


FAVO[®] S-TLIF CAGE FOR LUMBAR SPINAL FUSION

INSTRUMENTATION GUIDE



MADE IN GERMANY

TABLE OF CONTENTS

Preface	3
Indications / Contraindications	4
Access for the FAVO S-TLIF Cage	5
FAVO S-TLIF Cage – Instrumentation	7
Position correction of the FAVO S-TLIF Cage	17
PRODUCT INFORMATION	20
FAVO S-TLIF Implants	PI 02
FAVO S-TLIF Trial Implants	PI 05
FAVO S-TLIF and ROCCIA Instruments	PI 07
General instruments	PI 09
Alphabetical Index	PI 10

NOTE: This Guide describes the FAVO S-TLIF Cage instrumentation – it does not replace briefing by a physician experienced in the surgical technique of the spinal column. This guide should be treated as reference material rather than a prescribed surgical procedure.

We would be happy to assist you in finding a hospital that provides an opportunity to observe surgical procedures.


PREFACE

FAVO[®] S-TLIF CAGE – FOR LUMBAR SPINAL FUSION

The FAVO S-TLIF (straight transforaminal lumbar interbody fusion) Cage is an implant for primary stabilization and restoration of physiological curvature in the lumbar and thoracic spine. The cage is designed for transforaminal approaches.

The aim is to reduce discogenic back pain, correct deformities, remedy instabilities, restore intervertebral height, restore physiological lordosis and provide biomechanical support for bone fusion in the disc space.

It must be used with an additional posterior stabilization system. For posterior thoracic and lumbar stabilization, Silony recommends the use of a posterior spinal fixator (e.g., the VERTICALE system).



Indications

The FAVO S-TLIF Cage is indicated for intervertebral body fusion of the spine in skeletally mature patients. The system is designed for use with autogenous bone graft to facilitate fusion and supplemental internal spinal fixation systems cleared by the FDA for use in the thoracolumbar spine from T1-S1. The device is to be used in patients who have had at least six months of non-operative treatment.

The FAVO S-TLIF Cage is intended for use in interbody fusion at one or two contiguous levels in the thoracic spine from (T1-T12) and the thoracolumbar junction (T12-L1), following discectomy for the treatment of a symptomatic disc degeneration (DDD), including thoracic disc herniation (myelopathy and/or radiculopathy with or without axial pain).

The FAVO S-TLIF Cage is intended for use at one or two contiguous levels in the lumbar spine, from L1 to S1, for the treatment of degenerative disc disease (DDD) with up to Grade I spondylolisthesis. DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. The FAVO S-TLIF Cage can be used as an adjunct to fusion in patients diagnosed with multilevel degenerative scoliosis.

Contraindications

Contraindications include, but are not limited to:

- Anticipated or documented allergy or intolerance to the materials (e.g., titanium)
- Any case in which the chosen implants would be too large or too small to achieve a successful result
- Any patient for whom use of the implant would be in conflict with the anatomical structures
- Missing bone structures, that render good anchoring of the implants impossible (e.g., associated with fractures, tumors, osteoporosis, or infections)

NOTE: FAVO S-TLIF **must** be combined with additional stabilization. For posterior thoracic and lumbar intersegmental fusion procedures, Silony Medical recommends the use of a posterior spinal fixator (e.g. with the VERTICALE system).

CAUTION: Please carefully review the advice on indications and contraindications in the IFU of FAVO S-TLIF. The Instruction for use must also be carefully reviewed as it also contains important information that might lead to exclusion of the patient.

ACCESS FOR THE FAVO[®] S-TLIF CAGE

The access for the FAVO S-TLIF Cage is described below.

Position and approach

RI-1410*
ROCCIA Chisel 10 mm,
reinforced



The patient is positioned in the prone position as is common for the posterior approach. Exposing the abdomen decreases the load on the abdominal vessels. Corresponding bearing frames or padding underneath the pelvis and thorax can be used for this. The main incision is usually performed medially above the spinous processes depending on the spinal segments being treated. The spinal erector muscles are then displaced strictly subperiosteally on both sides and prepared until the anatomical structures of the spinal column are clearly exposed.

The transforaminal approach to the disc space is usually achieved by means of unilateral resection of the facet joint on the approach side (Fig. 1). The 6, 8 or 10 mm ROCCIA Chisels can be used for this, as well as standard instruments such as Luer forceps and swages.

* Representative of other chisels
see ROCCIA Instruments




Fig. 1 Resection of a facet joint with the chisel for a TLIF approach

NOTE: The FAVO S-TLIF System is compatible with the ROCCIA Instruments as described below.

FAVO[®] S-TLIF CAGE INSTRUMENTATION

Discectomy

RI-1020*
ROCCIA Ring Endplate Scraper
Straight



RI-1040**
ROCCIA Curette Straight



FI-4107
FAVO Shaver, 7 mm, parallel



GI-3101****
T-Handle



To begin with, the disc is incised with a standard scalpel. The disc material is loosened by means of shavers and is then removed using various standard forceps and available endplate scrapers and curettes (Fig. 2, Fig. 3). The fibrous ring is opened up in the process, and the nucleus and the inner fibrous ring are then removed and the surfaces of the end plates are roughened in order to prepare a spacious cage bed.

Various angled endplate scrapers are available to facilitate the removal of the intervertebral disc tissue in the wide lateral disc space.

- * Representative for other ring endplate scrapers (angled and curved)
- ** Representative for other curettes (curved)
- *** Representative for other shaver sizes
see ROCCIA and FAVO Instruments
- **** Representative for other T-handles
see General instruments

NOTE: We recommend to preserve the outer fibrous ring as support for the cage.




Fig. 2 Loosening of the disc material with a shaver




Fig. 3 Discectomy with curette

Preparing the disc space

RI-1020*
ROCCIA Ring Endplate Scraper
Straight



RI-1030
ROCCIA Box Endplate Scraper
Straight



RI-1050
ROCCIA Rasp Curved 45°



FI-4107**
FAVO Shaver, 7 mm, parallel



GI-3101***
T-Handle



For more extensive curettage, the ROCCIA Box Endplate Scraper is also available (Fig. 4). The surface of the remaining cartilaginous layer of the base plate and cover plate can be roughened with bone rasps, curettes and shavers (Fig. 5). The curved Endplate Scrapers in particular also facilitate the preparation of the opposite side.

- * Representative for other Ring Endplate Scrapers
- ** Representative for other Box Endplate Scrapers
- *** Representative for other Rasps
- **** Representative for other shaver sizes
see ROCCIA and FAVO Instruments

WARNING: FAVO Shavers are only suitable for mobilizing the disc material for subsequent removal and preparing the cover plates. They must not be used for distraction.

WARNING: Careful preparation of the disc space, especially extensive roughening of the end plates, provides the basis for better vascularization and successful bone fusion. Damage of the bony base and cover plate can lead to subsidence of the implant into the vertebral body.




Fig. 4 Box endplate scraper for extensive removal of disc material




Fig. 5 Rasp for roughening the cover plates

Distracting the disc space

FI-4207*
FAVO Paddle Sizer 7 mm,
parallel



GI-3101**
T-Handle



Blunt FAVO Paddle sizers with depth marking are available for distraction. The paddle sizers are connected to the T-handle by means of a quick-release system. For better orientation, the handle ends are aligned in the same way as the end of the paddle shavers. Two T-handles are available to enable rapid instrumentation.

For distraction, a blunt paddle sizer corresponding to the disc height is first inserted into the disc space and erected by 90° rotation (Fig. 6). The next distractors are then inserted in steps of 1 mm with the same movement until the desired height is achieved. The appropriate distraction height is achieved when the paddle sizer is under tension and provides a feeling of stability. Standard lamina spreader forceps can additionally be used for distraction.

* Representative for other paddle sizers
see FAVO Instruments

** Representative of other T-handles
see General instruments




Fig. 6 Spreading the disc space with the Paddle Sizer

CAUTION: Overdistraction should be avoided as it increases the risk of damaging the base plates and cover plates and subsequent subsidence of the implant, and jeopardizes the restoration of physiological curvature.

Selecting the trial sizer

FI-4207*
FAVO Paddle Sizer 7 mm,
parallel



GI-3101**
T-Handle



To determine the size of the disc space, the blunt paddle sizer with depth markings are available. Trial sizer can be selected on the basis of these measurements.

An appropriate trial sizer with 5° or 15° lordosis is available for each definitive cage size. The height of the trial sizer corresponds to the height of the final implant excluding tothing. Due to the tothing, the final implant is approx. 0.8 mm higher than the trial sizer.

* Representative for other paddle sizers,
see FAVO Instruments

** Representative for other T-handles
see General instruments

Determining the cage sizes with the trial sizers

FI-T07112805*
FAVO S-TLIF Trial Height 7mm,
5°



GI-3101**
T-Handle



The trial sizer combines all three possible implant lengths 28, 33 and 38 mm, in the respective height. The height of the implant should be similar to the normal disc height. The trial sizers are colourmarked analogously to the packaging of the implant to be implanted later on.

Should the image converter control require a larger lordosis, the sizer with 5° and 15° lordosis angle is also available.

Silony Medical recommends selecting the largest possible implant to maintain a maximal support surface including the anterior and posterior cortical area of the endplates. The trial sizer is carefully inserted into the intervertebral space by applying mild pressure (Fig. 7).

The length is determined in the image converter by means of X-ray markings on the trial sizer, which show the different lengths in the lateral X-ray image. (Fig. 8).

To ensure that the height of the intervertebral disc is preserved after loosening the distraction, the sizer must fit between the end plates after full distraction of the segment. Using the largest possible trial sizer for each individual patient maximizes the stability of the segment. If the trial sizer does not sufficiently fill in the intervertebral disc space, the next larger trial sizer must be used. If the trial sizer cannot be inserted because the intervertebral disc space is too small, either the next smaller size has to be used or the segment has to be distracted further using the aforementioned instruments.

* Representative for other trial sizes
see FAVO S-TLIF Trials




Fig. 7 Inserting the trial sizer




Fig. 8 Image converter control during insertion of the trial sizer


WARNING: Correctly selecting the cage size is critical to the success of the instrumentation and fusion.

Multitude of cage sizes

To optimize the treatment of the patient, a wide range of FAVO S-TLIF sizes is available.

The portfolio comprises sizes in 9 anterior heights (7 – 13 mm in 1-mm increments, and heights of 15mm and 17mm), three footprints (28x11 mm, 33x11 mm, 38x11 mm) and lordotic angles (5°, 10°, 15°) that allow the restoration of the sagittal profile.

The trial sizer with the particular colour code corresponds to the packaging of the implant to be implanted later on.



Filling of the cage

FI-4020
FAVO Inserter M4,
dismountable



FI-4050
FAVO S-TLIF Loading Block



FI-4051
FAVO S-TLIF Bone Graft Pusher



Remaining areas of the intervertebral disc space must be filled with autogenous bone graft either before or after implantation of the cage in order to achieve the largest possible fusion surface area.

The selected implant is screwed onto the respective FAVO Inserter, enabling the cage to be definitively inserted without requiring any further instrument change.




Fig. 9 Filling the cage with bone material in the loading block with pusher

CAUTION: The autogenous bone graft must be inserted well compressed into the cage.

WARNING: Proper filling of the disc space and implant are critical prerequisites for secure fusion. Use the provided loading block and bone graft pusher to ensure proper filling (Fig. 9).

Inserting the cage

FI-4020
FAVO Inserter M4,
dismountable



RI-1355
ROCCIA Slotted Mallet, solid



The trial sizer is removed just shortly before definitively placing the implant in order to prevent subsequent sinking.

For firm connection between cage and insertion instrument the turning handle of the insertion instrument is rotated in the direction of the closed lock (Fig. 10).

To release the cage, turn the knob in the direction of the open lock (Fig. 11).

The filled implant is carefully inserted into the disc space and the correct alignment of the implant is verified (Fig. 12). Slight pressure or careful hammering with the ROCCIA Slotted Mallet on the implant holder may be required.

Under X-ray control, the cage is inserted into the disc space. The neurogenic structures must be protected.




Fig. 10 Turning in the direction of the closed lock firmly connects the insertion instrument to the cage




Fig. 11 Turning in the direction of the open lock releases the cage from the insertion tool again




Fig. 12 Inserting the filled cage into the disc space

Optimal position of the FAVO® S-TLIF Cage

If possible, the inserter is left in the cage until an AP image and a lateral image with the image converter confirm the correct position of the cage.

The optimal position of the implanted FAVO S-TLIF cage (Fig. 13, Fig. 14) should be 35° +/-10 in the axial plane and positioned as anterior as is safe, close to the anterior vertebral endplate rim. In X-ray view, the optimal position and correct insertion angle can be confirmed by the implant located close to the anterior third of the disc space, the two lateral cage windows should display a clear shape and with the web in the middle.

The more anterior the cage is positioned, the better curvature can be achieved in the respective section of the spine. Avoid a position too far posterior. Bicortical support is ideal.

Once the cage has been successfully implanted, the remaining disc space must be filled up with autogenous bone graft to ensure secure fusion.




Fig. 13 Optimal position of the filled S-TLIF Cage




Fig. 14 X-ray control in the lateral projection

CAUTION: Due to the rough surface of the implant that provides good primary stability, the insertion in exact angle should be considered carefully to avoid corrections afterwards.

WARNING: The cage must be used with an additional posterior stabilization system (e.g., the VERTICALE system). The final steps of posterior fixation (e. g. insertion of the rod, compression and final tightening of the set screws) are completed after implantation of the cage.

NOTE: When implanting an S-TLIF Cage with 15° lordosis, it is imperative to provide the segment with additional stabilization.

CORRECTING THE POSITION OF THE FAVO S-TLIF CAGE

A straight implant driver and a hooked implant driver are provided for definitive positioning of the FAVO S-TLIF Cage. The use of these instruments is outlined below.

Correcting the position of the FAVO S-TLIF Cage

RI-1343
ROCCIA Hooked Implant Driver,
reinforced



RI-1355
ROCCIA Slotted Mallet, solid



RI-1340
ROCCIA Implant Driver Straight



FI-4020
FAVO Inserter M4,
dismountable



The small tip of the ROCCIA Hooked Implant Driver is carefully inserted through the posterior opening of the cage (Fig. 15). The cage can then be hit into the desired position with the slotted mallet. The hooked implant driver ensures a certain axial guidance during impaction.

The ROCCIA Implant Driver Straight is positioned straight onto the lateral posterior front of the cage (Fig. 16). The cage can then be carefully moved into the desired position with the slotted mallet.

The inserter has been designed to fulfil multiple functions. Using the rotating handle, the two instrument components can be disassembled so that both individual components are available for use. The main body can be attached to the cage after removal of the inner core and act as a driver by slightly hammering (Fig. 17). The cage can then be placed more anterior. The correcting should be performed under image converter control.




Fig. 15 Hooked Implant Driver for correcting the position of the cage




Fig. 16 Straight Implant Driver to finalize the position of the cage



Fig. 17 Main body of FAVO Inserter as a driver

Revision

FI-4020
FAVO Inserter M4,
dismountable



To remove the FAVO S-TLIF implant, the inner core of the FAVO insertion instrument can be used. The core is attached to the cage and tightened. The removal of a cage requires that the instrument be screwed in until it stops. The cage can be removed by slight careful hammering under the handle (Fig. 18).



Fig. 18. Revision of the FAVO S-TLIF Cage


FAVO[®] S-TLIF PRODUCT INFORMATION

FAVO [®] S-TLIF Implants by article number	PI 02
FAVO [®] S-TLIF Trial Implants by article number	PI 05
FAVO [®] S-TLIF and ROCCIA Instruments by article number	PI 06
General Instruments by article number	PI 09
FAVO [®] S-TLIF Alphabetical Index	PI 10

FAVO® S-TLIF Implants

Article number explanation for the cage, as an example

FAVO S-TLIF Cage, 13 x 11 x 28 mm, 5° lor.



Article number	Description	Illustration
S-FST-07112805-S	FAVO S-TLIF Cage 7 x 11 x 28 mm, 5° lor.	
S-FST-08112805-S	FAVO S-TLIF Cage 8 x 11 x 28 mm, 5° lor.	
S-FST-09112805-S	FAVO S-TLIF Cage 9 x 11 x 28 mm, 5° lor.	
S-FST-10112805-S	FAVO S-TLIF Cage 10 x 11 x 28 mm, 5° lor.	
S-FST-11112805-S	FAVO S-TLIF Cage 11 x 11 x 28 mm, 5° lor.	
S-FST-12112805-S	FAVO S-TLIF Cage 12 x 11 x 28 mm, 5° lor.	
S-FST-13112805-S	FAVO S-TLIF Cage 13 x 11 x 28 mm, 5° lor.	
S-FST-09112810-S	FAVO S-TLIF Cage 9 x 11 x 28 mm, 10° lor.	
S-FST-10112810-S	FAVO S-TLIF Cage 10 x 11 x 28 mm, 10° lor.	
S-FST-11112810-S	FAVO S-TLIF Cage 11 x 11 x 28 mm, 10° lor.	
S-FST-12112810-S	FAVO S-TLIF Cage 12 x 11 x 28 mm, 10° lor.	
S-FST-13112810-S	FAVO S-TLIF Cage 13 x 11 x 28 mm, 10° lor.	
S-FST-11281115-S	FAVO S-TLIF Cage 11 x 11 x 28 mm, 15° lor.	
S-FST-12281115-S	FAVO S-TLIF Cage 12 x 11 x 28 mm, 15° lor.	
S-FST-13281115-S	FAVO S-TLIF Cage 13 x 11 x 28 mm, 15° lor.	

System:
FAVO

Implant type:
S-TLIF

Configuration:
28 mm


Material:
Ti6Al4V ELI

FAVO® S-TLIF Implants

System:
FAVO
Implant type:
S-TLIF

Configuration:
33 mm

Material:
Ti6Al4V ELI

Article number	Description	Illustration
S-FST-08113305-S	FAVO S-TLIF Cage 8 x 11 x 33 mm, 5° lor.	
S-FST-09113305-S	FAVO S-TLIF Cage 9 x 11 x 33 mm, 5° lor.	
S-FST-10113305-S	FAVO S-TLIF Cage 10 x 11 x 33 mm, 5° lor.	
S-FST-11113305-S	FAVO S-TLIF Cage 11 x 11 x 33 mm, 5° lor.	
S-FST-12113305-S	FAVO S-TLIF Cage 12 x 11 x 33 mm, 5° lor.	
S-FST-13113305-S	FAVO S-TLIF Cage 13 x 11 x 33 mm, 5° lor.	
S-FST-15113305-S	FAVO S-TLIF Cage 15 x 11 x 33 mm, 5° lor.	
S-FST-10113310-S	FAVO S-TLIF Cage 10 x 11 x 33 mm, 10° lor.	
S-FST-11113310-S	FAVO S-TLIF Cage 11 x 11 x 33 mm, 10° lor.	
S-FST-12113310-S	FAVO S-TLIF Cage 12 x 11 x 33 mm, 10° lor.	
S-FST-13113310-S	FAVO S-TLIF Cage 13 x 11 x 33 mm, 10° lor.	
S-FST-15113310-S	FAVO S-TLIF Cage 15 x 11 x 33 mm, 10° lor.	
S-FST-12113315-S	FAVO S-TLIF Cage 12 x 11 x 33 mm, 15° lor.	
S-FST-13113315-S	FAVO S-TLIF Cage 13 x 11 x 33 mm, 15° lor.	
S-FST-15113315-S	FAVO S-TLIF Cage 15 x 11 x 33 mm, 15° lor.	

FAVO® S-TLIF Implants

Article number	Description	Illustration
S-FST-09113805-S	FAVO S-TLIF Cage 9 x 11 x 38 mm, 5° lor.	
S-FST-10113805-S	FAVO S-TLIF Cage 10 x 11 x 38 mm, 5° lor.	
S-FST-11113805-S	FAVO S-TLIF Cage 11 x 11 x 38 mm, 5° lor.	
S-FST-12113805-S	FAVO S-TLIF Cage 12 x 11 x 38 mm, 5° lor.	
S-FST-13113805-S	FAVO S-TLIF Cage 13 x 11 x 38 mm, 5° lor.	
S-FST-15113805-S	FAVO S-TLIF Cage 15 x 11 x 38 mm, 5° lor.	
S-FST-17113805-S	FAVO S-TLIF Cage 17 x 11 x 38 mm, 5° lor.	
S-FST-11113810-S	FAVO S-TLIF Cage 11 x 11 x 38 mm, 10° lor.	
S-FST-12113810-S	FAVO S-TLIF Cage 12 x 11 x 38 mm, 10° lor.	
S-FST-13113810-S	FAVO S-TLIF Cage 13 x 11 x 38 mm, 10° lor.	
S-FST-15113810-S	FAVO S-TLIF Cage 15 x 11 x 38 mm, 10° lor.	
S-FST-17113810-S	FAVO S-TLIF Cage 17 x 11 x 38 mm, 10° lor.	
S-FST-13113815-S	FAVO S-TLIF Cage 13 x 11 x 38 mm, 15° lor.	
S-FST-15113815-S	FAVO S-TLIF Cage 15 x 11 x 38 mm, 15° lor.	
S-FST-17113815-S	FAVO S-TLIF Cage 17 x 11 x 38 mm, 15° lor.	

System:
FAVO

Implant type:
S-TLIF

Configuration:
38 mm

Material:
Ti6Al4V ELI

FAVO® S-TLIF Trial Implants

System:
FAVO

Instrument type:
Trial implant

Configuration:
28-38 mm

Material:
Stainless Steel

Article number	Description	Illustration
FI-T07112805	FAVO S-TLIF Trial Height 7 mm, 5° lor.	
FI-T08112805	FAVO S-TLIF Trial Height 8 mm, 5° lor.	
FI-T09112805	FAVO S-TLIF Trial Height 9 mm, 5° lor.	
FI-T10112805	FAVO S-TLIF Trial Height 10 mm, 5° lor.	
FI-T11112805	FAVO S-TLIF Trial Height 11 mm, 5° lor.	
FI-T12112805	FAVO S-TLIF Trial Height 12 mm, 5° lor.	
FI-T13112805	FAVO S-TLIF Trial Height 13 mm, 5° lor.	
FI-T15112805	FAVO S-TLIF Trial Height 15 mm, 5° lor.	
FI-T17112805	FAVO S-TLIF Trial Height 17 mm, 5° lor.	

System:
FAVO

Instrument type:
Trial implant

Configuration:
28-38 mm


Material:
Stainless Steel

Article number	Description	Illustration
FI-T11112815	FAVO S-TLIF Trial Height 11 mm, 15° lor.	
FI-T12112815	FAVO S-TLIF Trial Height 12 mm, 15° lor.	
FI-T13112815	FAVO S-TLIF Trial Height 13 mm, 15° lor.	
FI-T15112815	FAVO S-TLIF Trial Height 15 mm, 15° lor.	
FI-T17112815	FAVO S-TLIF Trial Height 17 mm, 15° lor.	






FAVO® S-TLIF and ROCCIA® Instruments

Article number	Description	Illustration	Page
RI-1020	ROCCIA Ring Endplate Scraper Straight		8, 9
RI-1021	ROCCIA Ring Endplate Scraper Angled 25°		
RI-1022	ROCCIA Ring Endplate Scraper Curved 45°		
RI-1030	ROCCIA Box Endplate Scraper Straight		9
RI-1040	ROCCIA Curette Straight		8
RI-1041	ROCCIA Curette Curved Right 45°		8
RI-1042	ROCCIA Curette Curved Left 45°		8
RI-1050	ROCCIA Rasp Curved 45°		9
RI-1406	ROCCIA Chisel 6 mm Width, reinforced		6
RI-1408	ROCCIA Chisel 8 mm Width, reinforced		
RI-1410	ROCCIA Chisel 10 mm Width, reinforced		



FAVO® S-TLIF and ROCCIA® Instruments

Article number	Description	Illustration	Page
FI-4107	FAVO Shaver 7 mm, parallel		8, 9
FI-4108	FAVO Shaver 8 mm, parallel		
FI-4109	FAVO Shaver 9 mm, parallel		
FI-4110	FAVO Shaver 10 mm, parallel		
FI-4111	FAVO Shaver 11 mm, parallel		
FI-4112	FAVO Shaver 12 mm, parallel		
FI-4113	FAVO Shaver 13 mm, parallel		
FI-4114	FAVO Shaver 14 mm, parallel		
FI-4115	FAVO Shaver 15 mm, parallel		
FI-4116	FAVO Shaver 16 mm, parallel		
FI-4117	FAVO Shaver 17 mm, parallel		



FAVO® S-TLIF and ROCCIA® Instruments

Article number	Description	Illustration	Page		
FI-4207	FAVO Paddle Sizer 7 mm, parallel		10, 11		
FI-4208	FAVO Paddle Sizer 8 mm, parallel				
FI-4209	FAVO Paddle Sizer 9 mm, parallel				
FI-4210	FAVO Paddle Sizer 10 mm, parallel				
FI-4211	FAVO Paddle Sizer 11 mm, parallel				
FI-4212	FAVO Paddle Sizer 12 mm, parallel				
FI-4213	FAVO Paddle Sizer 13 mm, parallel				
FI-4214	FAVO Paddle Sizer 14 mm, parallel				
FI-4215	FAVO Paddle Sizer 15 mm, parallel				
FI-4216	FAVO Paddle Sizer 16 mm, parallel				
FI-4217	FAVO Paddle Sizer 17 mm, parallel		14, 15, 18, 19		
FI-4020	FAVO Inserter M4, dismountable				
RI-1330	ROCCIA Removal Adapter				no image
RI-1340	ROCCIA Implant Driver Straight				
RI-1343	ROCCIA Hooked Implant Driver, reinforced		18		
RI-1355	ROCCIA Slotted Mallet, solid				
			15, 18		

FAVO® S-TLIF Instruments

Article number	Description	Illustration	Page
FI-4051	FAVO S-TLIF Bone graft pusher		14
FI-4050	FAVO S-TLIF Loading Block		

General Instruments

Article number	Description	Illustration	Page
GI-2101	T-Handle, short		8, 9, 10, 11
GI-3101	T-Handle		

FAVO® S-TLIF Alphabetical Index

A-Z	Description	Article number	Page
B	Bone Graft Pusher	FI-4051	14
	Box Endplate Scraper Straight	RI-1030	9
C	Chisel 6 mm width, reinforced	RI-1406	6
	Chisel 8 mm width, reinforced	RI-1408	
	Chisel 10 mm width, reinforced	RI-1410	
	Curette, curved Left 45°	RI-1042	8
	Curette, curved Right 45°	RI-1041	
	Curette, straight	RI-1040	
H	Hooked Implant Driver, reinforced	RI-1343	18
I	Implant Driver Straight	RI-1340	18
	FAVO Inserter, M4, dismountable	FI-4020	14, 15, 18, 19
P	FAVO Paddle Sizer 7 mm, parallel	FI-4207	10, 11
	FAVO Paddle Sizer 8 mm, parallel	FI-4208	
	FAVO Paddle Sizer 9 mm, parallel	FI-4209	
	FAVO Paddle Sizer 10 mm, parallel	FI-4210	
	FAVO Paddle Sizer 11 mm, parallel	FI-4211	
	FAVO Paddle Sizer 12 mm, parallel	FI-4212	
	FAVO Paddle Sizer 13 mm, parallel	FI-4213	
	FAVO Paddle Sizer 14 mm, parallel	FI-4214	
	FAVO Paddle Sizer 15 mm, parallel	FI-4215	
	FAVO Paddle Sizer 16 mm, parallel	FI-4216	
	FAVO Paddle Sizer 17 mm, parallel	FI-4217	
R	Rasp, curved 45°	RI-1050	9
	Removal Adapter	RI-1330	No image
	Ring Endplate Scraper, curved 45°	RI-1022	8, 9
	Ring Endplate Scraper Straight	RI-1020	
	Ring Endplate Scraper, angled 25°	RI-1021	
S	Slotted Mallet, solid	RI-1355	15, 18
	FAVO Shaver 7 mm, parallel	FI-4107	8, 9
	FAVO Shaver 8 mm, parallel	FI-4108	
	FAVO Shaver 9 mm, parallel	FI-4109	
	FAVO Shaver 10 mm, parallel	FI-4110	
	FAVO Shaver 11 mm, parallel	FI-4111	
	FAVO Shaver 12 mm, parallel	FI-4112	
	FAVO Shaver 13 mm, parallel	FI-4113	
	FAVO Shaver 14 mm, parallel	FI-4114	
	FAVO Shaver 15 mm, parallel	FI-4115	
	FAVO Shaver 16 mm, parallel	FI-4116	
	FAVO Shaver 17 mm, parallel	FI-4117	

FAVO® S-TLIF Alphabetical Index

A-Z	Description	Article number	Page
T	T-Handle, short	GI-2101	8, 9, 10, 11
	T-Handle	GI-3101	
	FAVO S-TLIF Loading Block	FI-4050	14
	FAVO S-TLIF Trial Height 7 mm, 5° lor.	FI-T07112805	12
	FAVO S-TLIF Trial Height 8 mm, 5° lor.	FI-T08112805	
	FAVO S-TLIF Trial Height 9 mm, 5° lor.	FI-T09112805	
	FAVO S-TLIF Trial Height 10 mm, 5° lor.	FI-T10112805	
	FAVO S-TLIF Trial Height 11 mm, 5° lor.	FI-T11112805	
	FAVO S-TLIF Trial Height 12 mm, 5° lor.	FI-T12112805	
	FAVO S-TLIF Trial Height 13 mm, 5° lor.	FI-T13112805	
	FAVO S-TLIF Trial Height 15 mm, 5° lor.	FI-T15112805	
	FAVO S-TLIF Trial Height 17 mm, 5° lor.	FI-T17112805	
	FAVO S-TLIF Trial Height 11 mm, 15° lor.	FI-T11112815	
	FAVO S-TLIF Trial Height 12 mm, 15° lor.	FI-T1212815	
	FAVO S-TLIF Trial Height 13 mm, 15° lor.	FI-T13112815	
	FAVO S-TLIF Trial Height 15 mm, 15° lor.	FI-T15112815	
	FAVO S-TLIF Trial Height 17 mm, 15° lor.	FI-T17112815	



www.silony-medical.com

 **Silony Medical Corp.**
8200 NW 27th Street, Suite 104
Doral, FL 33122, USA
Phone +1 305 916 0016

 **Silony Medical GmbH**
Leinfelder Strasse 60
70771 Leinfelden-Echterdingen, Germany
Phone +49 711 78 25 25 0
Fax +49 711 78 25 25 11