


VERTICALE[®] CERVICAL SCREW ROD SYSTEM WITH OCCIPITO-CERVICAL FUSION

INSTRUMENTATION GUIDE



MADE IN GERMANY

TABLE OF CONTENTS

Preface	4
Position and approach	6
Opening the pedicle	6
Probing the pedicle	6
Drilling	7
Determining the screw dimensions	8
Tapping	9
Selection of pedicle screws	10
Preparing the pedicle screwdriver	11
Picking up the screws	12
Pedicle screw insertion	12
Countersinking the pedicle screw	13
Aligning the screw heads (tulip)	13
Selecting and sizing the rods	14
Inserting the rods	15
Temporarily tightening the set screw	15
Final tightening using the counter torque	16
Verification	16
Compression and distraction	17
Reduction with the reduction instrument	18
Fixing the rod in place with the reduction instrument	18
Size determination	19
Preparing the cross connectors	19
Inserting and temporarily fixing the cross connector	20
Final tightening using the cross connector counter torque	21
Preparing the lateral connector	22
Inserting and temporarily fixing the lateral connector	22
Final tightening using the counter torque	23
Determine plate position and shape	24
Drilling with occiput drills	25
Probing the occiput screw holes	26
Determining the occiput screw dimension	26
Tapping with the occiput tap	27
Selection of occiput screws	28
Inserting the occiput screws	28

TABLE OF CONTENTS

Selecting and sizing the rods for occiput	29
Inserting the rods	30
Temporarily tightening the set screw	30
Final tightening using the occiput counter torque	31
Verification	31
Connection for VERTICALE CERVICAL to VERTICALE (Thoracic Spine)	32
Inserting and tightening the VERTICALE CERVICAL rod connector inline	32
Inserting and tightening the VERTICALE CERVICAL rod connector triangle	33
Selecting and sizing the transition rods	34
Implant removal procedure	35
PRODUCT INFORMATION	
VERTICALE CERVICAL Implants	PI 03
VERTICALE CERVICAL Instruments	PI 10

NOTE: This guide describes the instrumentation of the VERTICALE Cervical posterior spinal fixation system. This guide does not replace briefing by a physician experienced in the instrumentation used in spinal surgery.
We would be happy to assist you in finding a hospital that provides an opportunity to observe surgical procedures.

PREFACE

VERTICALE[®] CERVICAL SCREW ROD SYSTEM INCLUDING OCCIPITO-CERVICAL FUSION

The VERTICALE CERVICAL System is a posterior double rod fixation system for immobilization and stabilization of the spinal segments of the craniocervical junction (occiput-C2), cervical spine (C1-C7), and thoracic spine (T1-T3). The system may only be used in the field of human medicine and consists of polyaxial screws, favored angle polyaxial (FA) screws, rods, occiput plates, connectors and the related instrumentation. Implants of the VERTICALE CERVICAL System are supplied sterile and are intended for single use.

Different screw placements will be possible:

- C1: Lateral mass screw
- C2: Isthmus screw, pedicle screw, laminar screw, transarticular
- C1-2: Length of current screw can be used for this connection
- C3-7: Lateral mass screw, pedicle screw, laminar screw
- T1-T3: Pedicle screw, laminar screw

NOTE: Ventral interbody support in the form of an intervertebral implant device, such as a Silony Cage System, is recommended for treating instabilities of the anterior spine and is used at the discretion of the operating surgeon and in accordance with the respective indication.

NOTE: Please also follow the Instructions for Use provided with each product. All instrumentation guides and Instructions for Use can be found on our eLabeling portal (<https://elabeling.silonymedical.com/>).

VERTICALE® CERVICAL INSTRUMENTATION GUIDE

In the following section, we begin by describing a monosegmental posterior VERTICALE CERVICAL Instrumentation with pedicle screws that forms the basis for all subsequent steps with additional instruments and implants. Multisegmental instrumentations or the implantation of other screws (e.g. lateral mass screws) are also performed according to these instructions.

Position and approach

The patient is positioned in the prone position as is common for the posterior approach. The skin incision is performed medially above the spinous processes corresponding to the spinal segments to be treated. The soft tissue is then dissected until the anatomical structures of the spinal column can be clearly seen.

Opening the pedicle

VI-0010*
VERTICALE CERVICAL Awl with
stop



VI-0020*
VERTICALE CERVICAL Probe,
straight



The desired screw insertion point into the pedicle is defined by means of anatomical landmarks and under X-ray control or other suitable inspection methods. The cortex is subsequently opened with the VERTICALE CERVICAL Awl with depth stop (Fig. 1). For safety reasons, the awl has a depth stop after 6 mm. The pedicle is opened further down to the cancellous bone of the vertebral body with the corresponding VERTICALE CERVICAL Probe.

* Additional instruments to open the pedicle can be found in the chapter "Instruments".




Fig. 1 Opening the pedicle with the awl or probe

Probing the pedicle

VI-0025
VERTICALE CERVICAL Pedicle
Feeler



The one-sided VERTICALE CERVICAL Pedicle Feeler can be used to check the prepared screw channel for possible perforations (Fig. 2).




Fig. 2 Probing the pedicle with the pedicle feeler

Drilling

VI-0230*
VERTICALE CERVICAL Drill Guide
Standard



VI-0256*
VERTICALE CERVICAL Drill Guide
XL



VI-0235*
VERTICALE CERVICAL Drill 2.4
mm standard



Various drill guides are available. A standard version for screw sizes between 10 mm and 30 mm, and two XL versions for screw lengths of 56 mm. Choose the corresponding drill guide for the desired screw dimension (determine the appropriate screw dimension by preoperative planning).

Insert the depth stop into the corresponding drill guide and set the depth stop to the required screw length. The adjustment of the length is made by pressing the button on the proximal part of the drill guide. The corresponding screw length is indicated by the laser marking on the drill guide (Fig. 3).

Separate drills are available for the respective drill guides.* The drills have diameters of 2.4 mm and 2.9 mm for 3.5 mm and 4.0 mm screws, respectively. The drills are color coded according to the anodization color of the screws, i.e. yellow for the 3.5 mm screws and blue for the 4.0 mm screws. Choose an appropriate drill for corresponding screw diameter ensuring the color coding of the drill matches the anodization colour of the screw.

The drills can be used with one of the VERTICALE CERVICAL quick-coupling handles or with a power tool.

Drill the hole into the pedicle or lateral mass to the appropriate depth and use intraoperative X-ray control if necessary (Fig. 4).

* Further drill options are shown in the chapter "Instruments".

NOTE: Confirm the correct adjustment of the depth stop and the desired drill length prior to drilling! Please correct the position of the depth stop if necessary.

NOTE: Ensure the sharpness of the drill prior to drilling! In case of blunt drill, use a new one.




Fig. 3 Assembling of the drill guide and insertion of the corresponding drill




Fig. 4 Drilling a hole into the pedicle

Determining the screw dimensions

VI-0030
VERTICALE CERVICAL Depth
Gauge



VI-0020*
VERTICALE CERVICAL Probe,
straight



Using the markings on the VERTICALE CERVICAL Probe, the dimension of the pedicle screw can be estimated. The markings are between 10 and 30 mm in increments of 10 mm. Use the VERTICALE CERVICAL Depth Gauge with markings between 8 and 56 mm with increments of 2 mm to confirm the required screw length (Fig. 5).

* Further probes are shown in the chapter "Instruments"





Fig. 5 Determining the length of the screws

NOTE: Do not bend the depth gauge or push by force into the bone as damage of the instrument may occur, which may lead to false length indications. Please check, whether the tip of the depth gauge is intact. Use the depth gauge only if the tip is intact and the laser marking is clearly visible.

Tapping

VI-0230*
VERTICALE CERVICAL Drill Guide
Standard



VI-0256
VERTICALE CERVICAL Drill Guide
XL



VI-0301**
VERTICALE CERVICAL T-Handle



VI-0035*
VERTICALE CERVICAL Tap
3.5 mm std



All VERTICALE CERVICAL screws have a self-tapping thread. However, for very hard bone structures (e.g., sclerotic bone) it may be necessary to pre-tap the thread into the bone. Taps for the screws with a diameter of 3.5 and 4.0 mm are available for this. All taps are color coded according to the anodization color of the screws, i.e., yellow for the 3.5 mm screws and blue for the 4.0 mm screws. Separate taps for the various drill guides are available. Choose a tap with suitable diameter and length ensuring the color coding of the tap matches the anodization color of the screw. The depth of tapping is controlled by the depth stop of the drill guides (Fig. 6). The taps can be used with one of the VERTICALE CERVICAL quick-coupling handles. After selecting the appropriate modular and cannulated handle (T-handle, long handle, with or without ratchet mechanism), it is connected to the corresponding VERTICALE CERVICAL Tap by locking into place (Fig. 7). The screw channel is prepared clockwise. The thread on the VERTICALE CERVICAL taps has a length of 16 mm. After cutting, the tap is disengaged by turning it counter clockwise. Cannulated taps are available for guided insertion using a guide-wire (Ø 1.3 mm). See appendix VERTICALE CERVICAL instruments.

* Further taps and drill guides are shown in the chapter "Instruments".

** Further handle options are shown in the chapter "Instruments".

NOTE: If using another guide wire than listed in the appendix, ensure that the length of the guide wire exceeds the length of implant, instrument and additional handle.

NOTE: Do not use a power tool for tapping!




Fig. 6 Inserting the depth stop and set the depth stop to the required length











Fig. 7 Tapping of screw holes

Selection of pedicle screws

To enable faster and easier identification, all VERTICALE CERVICAL screws are color coded by diameter. The lengths vary by 2 mm increments for the 3.5 mm and 4.0 mm screws. The 4.5 mm screw serves as a revision option and vary by 5 mm increments. Two different types of screw heads are available. The poly-head with an angulation of 40° in all directions and an FA-head (FA = far angle) with an increased angulation of 57° in either cranial or caudal direction. 4.0 mm and 4.5 mm screws are also available with cannulation. Moreover, a smooth shank option is available with 10 mm non-threaded screw shaft for the 4.0 mm screws.*

* Further screw options are shown in the chapter “Implants”.

	Ø 3.5 mm (10 - 30 mm à 2 mm)		Ø 4.0 mm (14 - 56 mm à 2 mm)		Ø 4.0 mm (20 - 36 mm à 2 mm)		Ø 4.5 mm (20 - 55 mm à 5 mm)	
10 mm								
12 mm								
14 mm								
16 mm								
18 mm								
20 mm								
22 mm								20 mm
24 mm								
26 mm								
28 mm								
30 mm								
32 mm								
34 mm								
36 mm								
38 mm								
40 mm								
42 mm								
44 mm								
46 mm								
48 mm								
50 mm								
52 mm								
54 mm								
56 mm								

NOTE: Using the A-P X-ray image, choose pedicle screws according to the pedicle diameter with the largest possible diameter. The length of the screw should be such that it reaches at least 2/3 of the diameter of the vertebral body, and in the best case the anterior edge of the vertebral body.

Preparing the pedicle screwdriver

VI-0130*
VERTICALE CERVICAL Pedicle SD
(screwdriver)



VI-0201**
VERTICALE CERVICAL Straight
handle, cannulated



The VERTICAL CERVICAL Pedicle Screwdriver is used to screw in the VERTICAL CERVICAL screws. It has to be assembled prior to use.

Mount the inner shaft of the screwdriver to the basic core by pushing the button at the proximal part of the basic core. To ensure better protection of the tissue, the screwdriver is equipped with a removable protection sleeve. It is attached, as shown, until it clicks into position (Fig. 8a).

The screwdriver is mounted onto the desired handle using the quick coupling on the handle. Different modular handles are available for use with the screwdrivers (with or without ratchet mechanism). Laser marking on the quick-coupling supports the right orientation of the modular handle for mounting.

* Further screwdrivers are shown in the chapter "Instruments".

** Further handle options are shown in the chapter "Instruments".



Fig. 8a Assembly of pedicle screwdriver

Picking up the screws

VI-0130*
VERTICALE CERVICAL Pedicle SD
(screwdriver)*



VI-0201**
VERTICALE CERVICAL Straight
handle, cannulated



Care must be taken to ensure an orthograde alignment between the tulip and screw shaft. The inner shaft of the VERTICALE CERVICAL Pedicle Screwdriver is first inserted deeply into the inner Torx of the screw shaft. After that, the threaded basic core is pushed towards the tulip and the internal thread of the tulip is connected to the external thread of the instrument by rotating to knob of the screwdriver clockwise and applying mild downward force with the instrument shaft (Fig. 9).

Ensure that the T-shaped tip is sufficiently inserted into the tulip for secure locking of the screw to the pedicle screwdriver.




Fig. 9 Attaching pedicle screwdriver to pedicle screw

Pedicle screw insertion

VI-0130*
VERTICALE CERVICAL Pedicle SD
(screwdriver)



VI-0201**
VERTICALE CERVICAL Straight
handle, cannulated



The VERTICALE CERVICAL pedicle screws are screwed into the prepared screw channel of the bone until the screw shaft is fully inserted into the pedicle (Fig. 10). Screwing too far into the pedicle can restrict the mobility of the tulip and make it difficult to insert the rod later. To disengage the instrument from the pedicle screw, maintain firm grip of screwdriver handle while rotating the rotation knob of sleeve counter-clockwise and pull instrument away from tulip. This process is repeated until all pedicle screws have been inserted. Verifying the correct positioning of the pedicle screws by means of an image intensifier in frontal and sagittal projection is strongly recommended.

* Further screwdrivers are shown in the chapter "Instruments"

** Further handle options are shown in the chapter "Instruments"




Fig. 10 Inserting the pedicle screw into the bone

NOTE: It is important that the polyaxiality of the tulip is not blocked. If necessary, the screw must be turned back a little.

Countersinking the pedicle screw

VI-0446
VERTICALE CERVICAL Pedicle
Screwdriver Ballhead



The VERTICAL CERVICAL Pedicle Screwdriver Ballhead is available for countersinking the pedicle screws (Fig. 11).




Fig. 11 Countersinking the pedicle screws

Aligning the screw heads (tulip)

VI-0350
VERTICALE CERVICAL Rod and
Tulip Adjuster



The VERTICAL CERVICAL screw heads are adjusted with the VERTICAL CERVICAL Rod and Tulip Adjuster. The adjuster is placed into the tulip and can then be used to align the tulip (by rotating and tilting) depending on how the rod will subsequently be inserted (Fig. 12).




Fig. 12 Aligning the pedicle screw heads with the rod and tulip adjuster

Selecting and sizing the rods

VI-0535
VERTICALE Phantom Rod 3.5 x
200mm



VI-0610
VERTICALE CERVICAL Rod
Bender



VI-0260
VERTICALE CERVICAL Rod
Cutter*



Various rod lengths with a diameter of 3.5 and 4.0 mm are available. Details can be found in chapter "Implants".

The VERTICALE CERVICAL Phantom Rod can be used to determine the required rod length and curvature. In order to estimate the required rod length, laser markings on the phantom rod are implemented in increments of 10 mm. Rods that are too long can be shortened with the VERTICALE CERVICAL Rod Cutter. Rotate the knurled wheel until the two arrows are aligned. Insert the rod into the corresponding hole (i.e. 3.5 mm or 4.0 mm). Repeatedly, squeeze the handles and initiate the ratcheting mechanism until the rod is cut. Before the next rod can be cut, rotate the knurled wheel until the two arrows are aligned again (Fig. 13a).

For individual anatomic adjustment of the rod, you can use the VERTICALE CERVICAL Rod Bender. Place the rod within the rod bender and squeeze the handles until the desired curvature is reached. Three different bending radii are possible (small, medium, large). In order to change the bending radius, pull the center knob and turn (~120°) to select the required bending radius. Make sure the center knob is arrested correctly and fixed in its position before contouring the rod (Fig. 13b).

* Further rod cutters are shown in the chapter "Instruments".

NOTE: The "cutting line" indicates the area where the rod will be cut. This line is located approximately 8 mm from the upper face of the instrument where the rod will be inserted.

NOTE: Any reverse bending of the rod decreases the integrity of the material and must be avoided. For this reason, bending of the rod should be performed gradually until the desired curvature is attained.




Fig. 13a Cutting the rod with the rod cutter



Fig. 13b Bending the rod with the rod bender

Inserting the rods

VI-0330
VERTICALE CERVICAL Rod
Holder



VI-0350
VERTICALE CERVICAL Rod and
Tulip Adjuster



The rods are inserted using the VERTICALE CERVICAL Rod Holder (Fig. 14).

If the rod is not placed deep enough into the pedicle tulip, it can be additionally maneuvered into the correct position with the VERTICALE Rod and Tulip Adjuster.

NOTE: The end of the rod shall visibly protrude from the last screw head to ensure sufficient contact between tulip, set screw and rod.




Fig. 14 Inserting the rods

Temporarily tightening the set screw

VI-0421
VERTICALE CERVICAL Setscrew
Starter double



The VERTICALE CERVICAL set screw is inserted with the double-sided VERTICALE CERVICAL Set Screw Starter. To do this, the Torx of the VERTICALE CERVICAL Set Screw Starter is equipped with a self-retaining geometry that holds the set screw in place during handling. Insert the set screw into the tulip of the pedicle screw. The rod is temporarily fixed by gently turning the set screw clockwise (Fig. 15).

NOTE: Set screws should always be inserted with a smooth clockwise rotation. To prevent tilting, a brief prior counter clockwise rotation can facilitate insertion of the set screw into the first thread.




Fig. 15 Inserting and temporarily tightening the set screw

Final tightening using the counter torque

VI-0440
VERTICALE CERVICAL Torque
Limiter 3 Nm



VI-0450
VERTICALE CERVICAL Counter
Torque



The VERTICAL CERVICAL Counter Torque is used to stabilize the rotation when tightening the VERTICAL CERVICAL Set Screw. In order to insert the set screw with guidance, the counter torque is placed directly onto the screw head. The VERTICAL CERVICAL Counter Torque can be conveniently mounted parallel or at right angles to the rod. The VERTICAL CERVICAL Torque Limiter (Torx 20) can then be guided by the counter torque and the set screw is tightened in its final position with a torque of 3 Nm (an audible click indicates that the torque has been reached). The same procedure must be repeated with all other set screws (Fig. 16). We recommend ensuring that the set screw is correctly seated by repeatedly tightening with the torque limiter. This is confirmed by two clicking sounds.




Fig. 16 Final tightening using the counter torque

Verification

Please check the correct position of the rods and screws as well as the result of compression and distraction maneuvers by final X-ray.

The result of the instrumentation is verified using images in two planes from an image intensifier.

VERTICALE® CERVICAL – INSTRUMENT-BASED REDUCTION AND CORRECTION OPTIONS

It is often necessary to perform intra-operative reductions and corrections of the implant devices with the help of VERTICALE CERVICAL Instruments in one or more segments.

Compression and distraction

VI-0620
VERTICALE CERVICAL
Distraction Pliers



VI-0630
VERTICALE CERVICAL
Compression Pliers



To compress or distract the pedicle screws, the VERTICALE CERVICAL Distraction or Compression Pliers are attached to the rod. The desired maneuver is performed by pressing together the respective pliers (Fig. 17 + 18). Both compression and distraction pliers feature a parallel mechanism for improved handling. The set screws are then final tightened with the VERTICALE CERVICAL Torque Limiter 3 Nm and the VERTICALE CERVICAL Counter Torque to secure the result of the compression or distraction maneuver.




Fig. 17 Compression with the compression pliers




Fig. 18 Distraction with the distraction pliers

NOTE: The set screws, at least one in the segment being corrected, must not be closed tightly during the maneuver.

Reduction with the reduction instrument

VI-0360
VERTICALE CERVICAL Reduction
Instrument



The VERTICALE CERVICAL Reduction Instrument is used to reduce the rod into the tulip of the pedicle screws. It is positively locked onto the designated hook groove (notch) at the verge of the tulip and the wings of the reduction instrument. The rod is then persuaded into the tulip by pressing together the handles of the reduction instrument. At the same time, the position of the vertebral body is corrected to posterior (Fig. 19).




Fig. 19 Reduction with the reduction instrument

NOTE: Consider both wings of the reduction instrument being in contact with the outer notch of the screw head (illustration) before pushing the rod into the tulip.

Fixing the rod in place with the reduction instrument

VI-0360
VERTICALE CERVICAL Reduction
Instrument



VI-0421
VERTICALE CERVICAL Setscrew
Starter double



Fixation of the rod is achieved using the VERTICALE CERVICAL Set Screw. It is inserted with the VERTICALE Set Screw Starter. To do this, the set screw is attached to the self-retaining Torx on the VERTICALE Set Screw Starter (Fig. 20). The rod is temporarily fixed by turning the set screw.

NOTE: Neither the set screw starter nor the reduction tool are designed for final screw tightening and may get damaged if applied torque is too high. For final tightening use the torque limiter and counter torque.




Fig. 20 Fixation with the reduction instrument and set screw starter

INSTRUMENTATION WITH THE VERTICALE[®] CERVICAL CROSS CONNECTOR

VERTICALE CERVICAL Cross Connectors (CC) are recommended to improve rotational stability, especially for instrumentations spanning multiple segments.

Size determination

VI-0830
VERTICALE CERVICAL Cross
Connector Sizer



Three different sizes of the VERTICAL CERVICAL Cross Connector implants are available with color-coding via the pre-mounted set screws.

In order to determine the required size, use the VERTICAL CERVICAL Cross Connector Indicator. For this, the instrument is attached to both rods (Fig. 21). Both rods have to lie completely in the two recesses of the instrument. Otherwise, another size range needs to be checked with the instrument.

NOTE: In case of two possible size indications please select the larger cross connector in order to prevent high strain and accidentally disassembly.




Fig. 21 Size determination of cross connector

Preparing the cross connectors

VI-0330
VERTICALE CERVICAL Rod
Holder



VI-0810
VERTICALE CERVICAL CC Torque
Limiter 2Nm



Ensure the pre-mounted set screws are loosened by using the VERTICAL CERVICAL CC Torque Limiter 2 Nm (Torx 15) prior to insertion onto the rods. Pay attention not to detach the set screws completely from the cross connector (Fig. 22).




Fig. 22 Preparing and inserting the cross connector

Inserting and temporarily fixing the cross connector

VI-0330
VERTICALE CERVICAL Rod
Holder



VI-0810
VERTICALE CERVICAL CC Torque
Limiter 2Nm



VI-0820
VERTICALE CERVICAL CC
Counter Torque



The cross connector can be engaged with the VERTICALE CERVICAL Rod Holder and inserted between the rods. Ensure that the cross connector is positioned correctly on the rod before tightening the set screws (Fig. 23). The cross connector can be adjusted in three dimensions, i.e. adjustment of the length, rotation around the transverse axis and rotation around the sagittal axis. After placing the cross connector on the rods, first the lateral set screws and then the medial set screw are temporarily fixed with the VERTICALE CERVICAL CC Torque Limiter 2 Nm by tightening the screws hand-tight.





Fig. 23 Inserting and temporarily fixing the cross connector

Final tightening using the cross connector counter torque

VI-0810
VERTICALE CERVICAL CC Torque
Limiter 2Nm



VI-0820
VERTICALE CERVICAL CC
Counter Torque



The VERTICALE CERVICAL CC Counter Torque is used to stabilize the rotation when tightening the cross connector set screw. The VERTICALE CERVICAL CC Torque Limiter (Torx 15) can then be guided by the counter torque and the cross connector set screws are tightened in their final position with a torque of 2 Nm (an audible click indicates that the torque has been reached). The same procedure must be repeated with all other cross connector set screws. Begin fixation with both lateral set screws prior to tightening the medial screw (Fig. 24).




Fig. 24 Final tightening with the counter torque

INSTRUMENTATION WITH THE VERTICALE[®] CERVICAL LATERAL CONNECTOR

Preparing the lateral connector

VI-0330
VERTICALE CERVICAL Rod
Holder



VI-0421
VERTICALE CERVICAL Setscrew
Starter double



After removing the lateral connector from the sterile packaging, the pre-assembled set screws must first be completely unscrewed counter-clockwise with the double-sided VERTICALE CERVICAL Set Screw Starter.

Inserting and temporarily fixing the lateral connector

VI-0330
VERTICALE CERVICAL Rod
Holder



VI-0421
VERTICALE CERVICAL Setscrew
Starter double



Place the lateral connector with the VERTICALE CERVICAL Rod Holder at the desired position between the VERTICALE CERVICAL rod and the VERTICALE CERVICAL pedicle screw. Use the VERTICALE CERVICAL Setscrew Starter double to first insert the set screw into the lateral connector and then into the tulip of the pedicle screw. The set screws are temporarily fixed with the VERTICALE CERVICAL Setscrew Starter double by tightening the screws hand-tight (Fig. 25).




Fig. 25 Pedicle screw with lateral connector

Final tightening using the counter torque

VI-0440
VERTICALE CERVICAL Torque
Limiter 3 Nm



VI-0450
VERTICALE CERVICAL Counter
Torque



The VERTICALE CERVICAL Counter Torque is used to stabilize the rotation when tightening the VERTICALE CERVICAL Set Screw. In order to insert the set screw with guidance, the counter torque is placed directly onto the screw head. The VERTICALE CERVICAL Counter Torque can be conveniently mounted parallel or at right angles to the rod. The VERTICALE CERVICAL Torque Limiter (Torx 20) can then be guided by the counter torque and the set screw is tightened in its final position with a torque of 3 Nm (an audible click indicates that the torque has been reached). The same procedure must be repeated with all other set screws.

We recommend ensuring that the set screw is correctly seated by repeatedly tightening with the torque limiter. This is confirmed by two clicking sounds.

VERTICALE® CERVICAL OCCIPUT INSTRUMENTATION

Determine plate position and shape

VI-0700
VERTICALE CERVICAL Occiput
Plate Holder



VI-0710
VERTICALE CERVICAL OC Plate
Bender right



VI-0711
VERTICALE CERVICAL OC Plate
Bender left



The VERTICALE CERVICAL Occiput Plate Holder is used to determine the appropriate plate positioning and screw placement. Three different designs of occiput plates are available. Choose the appropriate plate design for the patient's anatomy.

If required, use the plate bender to bend the plate to the desired shape. For bending the lateral wings of the plate insert them into the proximal part of the plate bender (Fig. 26). A self-retaining mechanism prevents the plate from falling down during bending. The area of the anchors can be bent by placing the distal flat end of the plate bender around the anchors. Gently bend the plate to the desired shape. The bending should be performed only in the bending zones to avoid damaging the area of the screw holes.

NOTE: The appropriate plate size and plate position should be determined by preoperative imaging techniques.

NOTE: Do not bend the plate more than 15°. Do not perform reverse bending of the plate in order not to decrease the fatigue life of the material.





Fig. 26 Bending of the desired plate shape

Drilling with occiput drills

VI-0700
VERTICALE CERVICAL Occiput
Plate Holder



VI-0720
VERTICALE CERVICAL OC Drill
Guide



VI-0731*
VERTICALE CERVICAL OC Drill,
flexible



Insert the depth stop into the drill guide and set the depth stop to the required screw length. The adjustment of the length is made by pressing the button on the proximal part of the drill guide (Fig. 27). The corresponding screw length is indicated by the laser marking on the drill guide (Determine the appropriate screw dimension by preoperative planning). Drills are available in both straight and flexible versions (with a universal joint for difficult anatomy).^{*} All drills have a diameter of 3.4 mm. The drills can be used with one of the VERTICALE CERVICAL quick-coupling handles or with a power tool. If a power drill is used in combination with the flexible drill, avoid angulation of more than 45° of the universal joint during drilling as damage of the drill may occur.

Place the tip of the drill guide in one of the holes for the occiput plate. Drill the initial occipital pilot hole in one of the midline holes to the appropriate depth and use intraoperative X-ray control if necessary (Fig. 28).

^{*} Further drill options are shown in the chapter "Instruments".

NOTE: Confirm the correct adjustment of the depth stop and the desired drill length prior to drilling! Please correct the position of the depth stop if necessary.

NOTE: Ensure the sharpness of the drill prior to drilling! In case of blunt drill, use a new one.




Fig. 27 Assemble the Drill Guide




Fig. 28 Drilling the occiput screw holes

Probing the occiput screw holes

VI-0025
VERTICALE CERVICAL Pedicle
Feeler



The one-sided VERTICAL CERVICAL Pedicle Feeler can be used to check the prepared screw channel for possible perforations (Fig. 29).




Fig. 29 Probing the occipital screw holes

Determining the occiput screw dimension

VI-0030
VERTICALE CERVICAL Depth
Gauge



Use the VERTICAL CERVICAL Depth Gauge with markings between 8 and 56 mm with increments of 2 mm to confirm the required screw dimension.

NOTE: Do not bend the depth gauge or push by force into the bone as damage of the instrument may occur, which may lead to false length indications. Please check, whether the tip of the depth gauge is intact. Use the depth gauge only if the tip is intact and the laser marking is clearly visible.

Tapping with the occiput tap

VI-0700
VERTICALE CERVICAL Occiput
Plate Holder



VI-0720
VERTICALE CERVICAL OC Drill
Guide



VI-0741*
VERTICALE CERVICAL OC Taps,
flexible



VI-0301**
VERTICALE CERVICAL T-Handle



For very hard bone structures (e.g. sclerotic bone) it may be necessary to pre-tap the thread into the bone. Taps are available for screws with a diameter of 4.5 mm. The depth of the tapping is controlled by the depth stop of the drill guides (compare to chapter "Drilling with occiput drills"). The taps can be used with one of the VERTICALE CERVICAL quick-coupling handles. After selecting the appropriate handle, it is connected to the corresponding VERTICALE CERVICAL Tap by locking into place (Fig. 30). The screw channel is prepared clockwise. The thread on the VERTICALE CERVICAL Tap has a length of 16 mm. After cutting, the tap is disengaged by turning it counter clockwise.

* Further taps are shown in the chapter "Instruments".

** Further handle options are shown in the chapter "Instruments".




Fig. 30 Tapping the occiput screw holes

NOTE: Do not use power tool for tapping!

Selection of occiput screws

To enable faster and easier identification, all VERTICALE CERVICAL occiput screws are color coded by diameter, i.e., violet and silver for the 4.5 and 5.2 mm screws, respectively. The lengths vary between 4 and 20 mm by 2 mm increments. The 5.2 mm screw serves as a revision option.

Inserting the occiput screws

VI-0700
VERTICALE CERVICAL Occiput
Plate Holder



VI-0751*
VERTICALE CERVICAL OC Screw-
driver, flexible



VI-0761
VERTICALE CERVICAL OC Count-
er Torque



The VERTICALE CERVICAL occiput bone screw is inserted with the VERTICALE CERVICAL Occiput Screwdriver. To do this, the torx of the screwdriver is equipped with a self-retaining geometry that holds the occiput screw in place during handling. Insert the occiput bone screw into the prepared bone channel. The occiput bone screw is temporarily fixed by gently turning the screwdriver clockwise. For difficult anatomy, a flexible screwdriver with universal joint is also available. In order to support the guidance of the flexible screwdriver, the occiput counter torque can be used as a guide (Fig. 31).

Insert all other occiput bone screws in the same manner. Then, final tightening is performed with one of the occiput screwdrivers.

* Further occiput screwdriver options are shown in the chapter "Instruments".





Fig. 31 Inserting the occiput screws

Selecting and sizing the rods for occiput

VI-0535
VERTICALE CERVICAL Phantom
Rod



VI-0610
VERTICALE CERVICAL Rod
Bender



VI-0611
VERTICALE CERVICAL Rod
Bending Tube



VI-0260*
VERTICALE CERVICAL Rod Cutter



Despite the standard rods, prebent rods are available especially for the occiput region with an angulation of 55° and two different diameters, i.e., 3.5 and 4.0 mm. The length of the prebent rod is 200 mm for the cervical part and 90 mm for the occipital part (Details can be found in the chapter “Implants”). For further individual anatomic adjustment of the rod, the VERTICALE CERVICAL Rod Bender or the VERTICALE CERVICAL Rod Bending Tubes can be used. Rods that are too long can be shortened with the VERTICALE CERVICAL Rod Cutter.

The VERTICALE CERVICAL Phantom Rod can be used to determine the required rod length and curvature. In order to estimate the required rod length, laser markings on the phantom rod are implemented in increments of 10 mm.

* Further rod cutters are shown in the chapter “Instruments”.

NOTE: Any reverse bending of the rod decreases the integrity of the material and must be avoided. For this reason, bending of the rod should be performed gradually until the desired curvature is attained.

NOTE: When using the VERTICALE CERVICAL bending tubes, do not make sharp bends or reverse bending of the rods in order to avoid potential fatigue of the implant.

Inserting the rods

VI-0330
VERTICALE CERVICAL Rod
Holder



The VERTICALE CERVICAL Rod Holder is used to place the bent rod into the anchors of the occiput plate as well as the tulips of the pedicle screws (Fig. 32).




Fig. 32 Inserting the rods

Temporarily tightening the set screw

VI-0750
VERTICALE CERVICAL OC
Screwdriver



The VERTICALE CERVICAL Set Screw is inserted with the VERTICALE CERVICAL Screwdriver. Insert the set screw into the anchor of the occiput plate. The rod is temporarily fixed by gently turning the set screw clockwise (Fig. 33).




Fig. 33 Tightening the set screws

NOTE: Set screws should always be inserted with a smooth clockwise rotation. To prevent tilting, a brief prior counter clockwise rotation can facilitate insertion of the set screw into the first thread.

Final tightening using the occiput counter torque

VI-0760
VERTICALE CERVICAL OC Torque
Limiter 3 Nm



VI-0761
VERTICALE CERVICAL OC
Counter Torque



The VERTICALE CERVICAL Occiput Counter Torque is used to stabilize the rotation when tightening the VERTICALE CERVICAL set screw. In order to insert the set screw with guidance, the counter torque is placed directly onto the anchor of the occiput plate. The VERTICALE CERVICAL Occiput Counter Torque can be comfortably mounted parallel or at right angles to the rod. The VERTICALE CERVICAL Occiput Torque Limiter (Torx 20) can then be guided by the counter torque and the set screw is tightened in its final position with a torque of 3 Nm (an audible click indicates that the torque has been reached). The same procedure must be repeated with the other set screw within the anchor of the occiput plate (Fig. 34). We recommend ensuring that the set screw is correctly seated by repeatedly tightening with the torque limiter. This is confirmed by two clicking sounds.




Fig. 34 Final tightening of the set screws

Verification

Please check the correct position of the plate, rods and screws by final X-ray.

The result of the instrumentation is verified using images in two planes from an image intensifier.

CONNECTION OF VERTICALE® CERVICAL TO VERTICALE (THORACIC SPINE)

In order to combine the VERTICALE CERVICAL System with the VERTICALE System in the thoracic spine two options are available: Two different types of inline connectors as well as transition rods with different dual-diameters.

Inserting and tightening the rod connector inline

VI-1830
VERTICALE Rod and Cross
Connector Holder*



VI-1810
VERTICALE T20 Screwdriver
7 Nm*



VERTICALE CERVICAL Rod Connectors inline are available for connecting the VERTICALE CERVICAL 3.5 or 4.0 mm rods lengthwise in order to enable a connection to be made between a new segment and a previously treated fusion segment using 5.5 mm rods.

The VERTICALE Rod and Cross Connector Holder is used to engage the VERTICALE CERVICAL Rod Connector Inline. Before the set screws are tightened, the 5.5 mm rod has to be inserted into the larger opening of the VERTICALE CERVICAL Rod Connector Inline as deeply as possible. This can be checked in the viewing panel between the two set screws. Then, either the 3.5 or the 4.0 mm VERTICALE CERVICAL rod is inserted into the smaller opening of the VERTICALE CERVICAL Rod Connector Inline. The final fixing into place is done by tightening the set screws with the VERTICALE T20 Screwdriver 7 Nm to a specified torque of 7 Nm (an audible click indicates that the torque has been reached). The same procedure must be repeated with all other set screws.

We recommend ensuring that the screw is correctly seated by repeatedly tightening with the torque limiter. This is confirmed by two clicking sounds.

* Instrument contained in VERTICALE Open Trays

Inserting and tightening the rod connector triangle

VI-1820*
VERTICALE Domino Holder



VI-1810*
VERTICALE T20 Screwdriver
7 Nm



The VERTICALE Domino Holder is used to engage the VERTICALE CERVICAL Triangle Rod Connector. Before the set screws are tightened, the 5.5 mm rod has to be inserted into the larger opening of the VERTICALE CERVICAL Rod Connector Triangle until both set screws have contact to the rod. Then, either the 3.5 or the 4.0 mm rod is inserted into the smaller opening of the VERTICALE CERVICAL Rod Connector Triangle. The final fixing into place is done by tightening the set screws with the VERTICALE T20 Screwdriver 7 Nm to a specified torque of 7 Nm (an audible click indicates that the torque has been reached). The same procedure must be repeated with all other set screws.

We recommend ensuring that the screw is correctly seated by repeatedly tightening with the torque limiter. This is confirmed by two clicking sounds.

* Instrument contained in VERTICALE Open Trays

NOTE: For using the VERTICALE CERVICAL Rod Connector (Inline and Triangle), ensure that the end of the overlapping rod is long enough to be inserted into the connectors.

NOTE: For the connection of the VERTICALE CERVICAL Rod Connectors the VERTICALE Open Trays are required and the final tightening has to be performed according to the VERTICALE Open Instrumentation Guide.

Selecting and sizing the transition rods

VI-0535
VERTICALE CERVICAL Phantom
Rod 3.5 x 200 mm



VI-0610
VERTICALE CERVICAL Rod
Bender



VI-0260*
VERTICALE CERVICAL Rod Cutter



VI-1261**
VERTICALE Rod Cutter, 5.5 mm



VI-1270**
VERTICALE French Bender



The VERTICALE CERVICAL System offers two dual diameter (transition) rods that can be connected to a thoracolumbar rod screw system with 5.5 mm rods, i.e., 3.5 and 5.5 mm or 4.0 and 5.5 mm. Details can be found in the chapter “Implants”.

Individual anatomic adjustment of the 3.5 and 4.0 mm diameter part of the transition rod can be performed according to the chapter: “Selecting and sizing the rods”.

* Further rod cutters are shown in the chapter “Instruments”.

** Instrument contained in VERTICALE Open Trays

NOTE: For the rod bending and rod cutting steps for the 5.5 mm rods the VERTICALE Open Tray is required and the final tightening of the 5.5 mm rods has to be performed according to the VERTICALE Open Instrumentation Guide.

NOTE: Any reverse bending of the rod decreases the integrity of the material and must be avoided. For this reason, bending of the rod should be performed gradually until the desired curvature is attained.

Implant removal procedure

To remove an implant, please execute the following steps as described. Pay special attention during the working steps on the loosened implants and screws.

Step 1: Removal of the Cross Connector

Use the VERTICALE CERVICAL CC Torque Limiter 2 Nm (VI-0810) to loosen the set screw on both sides of the cross connector and the central position by turning the handle counterclockwise. The VERTICALE CERVICAL CC Counter Torque (VI-0820) should be used to stabilize the rotation when loosening the set screws. Once the screws are loosened, use the VERTICALE CERVICAL Rod Holder (VI-0330) to remove the cross connector from the rods.

Step 2: Removal of the Set Screws

Insert the VERTICALE CERVICAL Torque Limiter 3 Nm (VI-0440) and turn the set screw counterclockwise until it loosens. The VERTICALE CERVICAL Counter Torque (VI-0450) should be used to stabilize the rotation when loosening the set screw. All VERTICALE CERVICAL Torque Limiters are equipped with a self-retaining mechanism at their tip. Therefore, the set screws are attached to the instrument and can be securely removed from the pedicle screws. Repeat this step until all set screws have been removed.

Step 3: Removal of the Lateral Offset Connector

To remove a lateral offset connector ensure that the set screw of the connected pedicle or lateral-mass screw has already been removed. Otherwise remove the set screw as described in step 2.

To remove the lateral offset connector (tulip design), remove the set screw as described in Step 2. Once the set screw is removed, use the VERTICALE CERVICAL Rod Holder (VI-0330) to remove the lateral offset connector from the rod and the pedicle screw.

Step 4: Removal of the Rods/ Inline Connector

Once all of the set screws have been removed, take the rod with the VERTICALE CERVICAL Rod Holder (VI-0330) and lift it up to remove the rod from the screw heads.

To remove a rod-to-rod connector (Inline, Triangle), use the VERTICALE T20 Screwdriver (VI-1810) or the VERTICAL Screwdriver Shaft (VI-1446) to loosen the set screws by turning them counterclockwise.

Use the VERTICALE Bar and Cross Connector Forceps (VI-1830) or the VERTICALE Domino Forceps (VI-1820) to remove the rod-to-rod connector.

Step 5: Removal of the Pedicle Screws

The VERTICALE CERVICAL Pedicle SD X15 (VI-0130) or the VERTICALE CERVICAL PSD Inner Shaft (VI-0130.3) is used to remove the VERTICALE CERVICAL pedicle screws. For attachment of the pedicle screw, the VERTICALE CERVICAL Pedicle Screwdriver is fully inserted into the inner Torx of the screw shaft and rotated inwards via the sleeve of the pedicle screwdriver. Rotate the screwdriver counterclockwise until the screw is fully out of the bone.

You can also use the VERTICALE CERVICAL Pedicle SD Ballhead (VI-0446) to remove the pedicle screws.

Step 6: Removal of the Occiput Plate

Use the VERTICALE CERVICAL OC Screwdriver (VI-0750) or the VERTICALE CERVICAL OC Screwdriver, flex (VI-0751) to remove the set screws from the plate tulips. Remove the rod as described in step 4.

Use the VERTICALE CERVICAL OC Screwdriver (VI-0750) or the VERTICALE CERVICAL OC Screwdriver, flex (VI-0751) also for removing the occiput screws. Before removing the last occiput screw take the VERTICALE CERVICAL Occiput Plate Holder (VI-0700) to hold the plate to prevent from traveling.


VERTICALE® CERVICAL PRODUCT INFORMATION

VERTICALE CERVICAL Implants by article number	PI 02 – 13
VERTICALE CERVICAL Instruments by article number	PI 10 – 13

VERTICALE[®] CERVICAL Implants


Article number explanation for screws, as examples

VERTICALE CERVICAL Poly Screw Ø 3.5 x 20 mm, solid




Diameter


Differentiation by color coding




Shaft dimensions




Implant type – Axial Mobility



Shaft configuration – Shape



VERTICALE CERVICAL 3 Occi Screw 4.5 x 8 mm




Diameter

Differentiation by color coding



Shaft dimensions



VERTICALE® CERVICAL Occiput Implants

System:
VERTICALE CERVICAL

Implant type:
Occiput screw

Material:
Ti6Al4V ELI

All articles are sterile
packed.


- 2:2 pieces in one package

- 3:3 pieces in one package



Article number	Description	Figure
S-VOS-4504-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 4 mm	
S-VOS-4506-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 6 mm	
S-VOS-4508-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 8 mm	
S-VOS-4510-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 10 mm	
S-VOS-4512-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 12 mm	
S-VOS-4514-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 14 mm	
S-VOS-4516-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 16 mm	
S-VOS-4518-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 18 mm	
S-VOS-4520-3	VERTICALE CERVICAL Occi 3 Screw 4.5 x 20 mm	
S-VOS-5204-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 4 mm	
S-VOS-5206-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 6 mm	
S-VOS-5208-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 8 mm	
S-VOS-5210-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 10 mm	
S-VOS-5212-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 12 mm	
S-VOS-5214-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 14 mm	
S-VOS-5216-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 16 mm	
S-VOS-5218-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 18 mm	
S-VOS-5220-3	VERTICALE CERVICAL Occi 3 Screw 5.2 x 20 mm	
S-VOS-4504-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 4 mm	
S-VOS-4506-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 6 mm	
S-VOS-4508-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 8 mm	
S-VOS-4510-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 10 mm	
S-VOS-4512-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 12 mm	
S-VOS-4514-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 14 mm	
S-VOS-4516-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 16 mm	
S-VOS-4518-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 18 mm	
S-VOS-4520-2	VERTICALE CERVICAL Occi 2 Screw 4.5 x 20 mm	

VERCALE® CERVICAL Occiput Implants

Article number	Description	Figure
S-VOS-5204-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 4 mm	
S-VOS-5206-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 6 mm	
S-VOS-5208-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 8 mm	
S-VOS-5210-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 10 mm	
S-VOS-5212-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 12 mm	
S-VOS-5214-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 14 mm	
S-VOS-5216-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 16 mm	
S-VOS-5218-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 18 mm	
S-VOS-5220-2	VERCALE CERVICAL Occi 2 Screw 5.2 x 20 mm	

System:
VERCALE CERVICAL

Implant type:
Occiput screw




Material:
Ti6Al4V ELI

All articles are sterile packed.

- 2:2 pieces in one package

- 3:3 pieces in one package

⚠ D60035

Article number	Description	Figure
S-VOC-LAT-5	VERCALE CERVICAL Occipital Pl. Lat. 5	
S-VOC-LAT-4	VERCALE CERVICAL Occipital Pl. Lat. 4	
S-VOC-MID-3	VERCALE CERVICAL Occipital Pl. Mid. 3	

System:
VERCALE CERVICAL

Implant type:
Occiput plate

Material:
Ti6Al4V ELI

All articles are sterile packed and include two set screws for rod attachment.

⚠ D60035


VERTICALE® CERVICAL Implants

System:
VERTICALE CERVICAL


Implant type:
Polyaxial screw

Material:
Ti6Al4V ELI

All articles are sterile packed
and include a set screw.



 D60035

Article number	Description	Figure
S-VCP-3510-S	VERTICALE CERVICAL Poly Scr 3.5 x 10 mm sol	
S-VCP-3512-S	VERTICALE CERVICAL Poly Scr 3.5 x 12 mm sol	
S-VCP-3514-S	VERTICALE CERVICAL Poly Scr 3.5 x 14 mm sol	
S-VCP-3516-S	VERTICALE CERVICAL Poly Scr 3.5 x 16 mm sol	
S-VCP-3518-S	VERTICALE CERVICAL Poly Scr 3.5 x 18 mm sol	
S-VCP-3520-S	VERTICALE CERVICAL Poly Scr 3.5 x 20 mm sol	
S-VCP-3522-S	VERTICALE CERVICAL Poly Scr 3.5 x 22 mm sol	
S-VCP-3524-S	VERTICALE CERVICAL Poly Scr 3.5 x 24 mm sol	
S-VCP-3526-S	VERTICALE CERVICAL Poly Scr 3.5 x 26 mm sol	
S-VCP-3528-S	VERTICALE CERVICAL Poly Scr 3.5 x 28 mm sol	
S-VCP-3530-S	VERTICALE CERVICAL Poly Scr 3.5 x 30 mm sol	

VERTICALE® CERVICAL Implants

Article number	Description	Figure
S-VCP-4014-K	VERTICALE CERVICAL Poly Scr 4.0 x 14 mm can	
S-VCP-4016-K	VERTICALE CERVICAL Poly Scr 4.0 x 16 mm can	
S-VCP-4018-K	VERTICALE CERVICAL Poly Scr 4.0 x 18 mm can	
S-VCP-4020-K	VERTICALE CERVICAL Poly Scr 4.0 x 20 mm can	
S-VCP-4022-K	VERTICALE CERVICAL Poly Scr 4.0 x 22 mm can	
S-VCP-4024-K	VERTICALE CERVICAL Poly Scr 4.0 x 24 mm can	
S-VCP-4026-K	VERTICALE CERVICAL Poly Scr 4.0 x 26 mm can	
S-VCP-4028-K	VERTICALE CERVICAL Poly Scr 4.0 x 28 mm can	
S-VCP-4030-K	VERTICALE CERVICAL Poly Scr 4.0 x 30 mm can	
S-VCP-4032-K	VERTICALE CERVICAL Poly Scr 4.0 x 32 mm can	
S-VCP-4034-K	VERTICALE CERVICAL Poly Scr 4.0 x 34 mm can	
S-VCP-4036-K	VERTICALE CERVICAL Poly Scr 4.0 x 36 mm can	
S-VCP-4038-K	VERTICALE CERVICAL Poly Scr 4.0 x 38 mm can	
S-VCP-4040-K	VERTICALE CERVICAL Poly Scr 4.0 x 40 mm can	
S-VCP-4042-K	VERTICALE CERVICAL Poly Scr 4.0 x 42 mm can	
S-VCP-4044-K	VERTICALE CERVICAL Poly Scr 4.0 x 44 mm can	
S-VCP-4046-K	VERTICALE CERVICAL Poly Scr 4.0 x 46 mm can	
S-VCP-4048-K	VERTICALE CERVICAL Poly Scr 4.0 x 48 mm can	
S-VCP-4050-K	VERTICALE CERVICAL Poly Scr 4.0 x 50 mm can	
S-VCP-4052-K	VERTICALE CERVICAL Poly Scr 4.0 x 52 mm can	
S-VCP-4054-K	VERTICALE CERVICAL Poly Scr 4.0 x 54 mm can	
S-VCP-4056-K	VERTICALE CERVICAL Poly Scr 4.0 x 56 mm can	
S-VCP-4520-K	VERTICALE CERVICAL Poly Scr 4.5 x 20 mm can	
S-VCP-4525-K	VERTICALE CERVICAL Poly Scr 4.5 x 25 mm can	
S-VCP-4530-K	VERTICALE CERVICAL Poly Scr 4.5 x 30 mm can	
S-VCP-4535-K	VERTICALE CERVICAL Poly Scr 4.5 x 35 mm can	
S-VCP-4540-K	VERTICALE CERVICAL Poly Scr 4.5 x 40 mm can	
S-VCP-4545-K	VERTICALE CERVICAL Poly Scr 4.5 x 45 mm can	
S-VCP-4550-K	VERTICALE CERVICAL Poly Scr 4.5 x 50 mm can	
S-VCP-4555-K	VERTICALE CERVICAL Poly Scr 4.5 x 55 mm can	


System:
VERTICALE CERVICAL

Implant type:
Polyaxial screw

Typing:
cannulated

Material:
Ti6Al4V ELI

All articles are sterile
packed and include a set
screw.



 D60035


VERTICALE® CERVICAL Implants

System:
VERTICALE CERVICAL

Implant type:
Far angle screw
(FA screw)

Material:
Ti6Al4V ELI

All articles are sterile
packed and include a set
screw.



Article number	Description	Figure
S-VCF-3510-S	VERTICALE CERVICAL FA Scr 3.5 x 10 mm solid	
S-VCF-3512-S	VERTICALE CERVICAL FA Scr 3.5 x 12 mm solid	
S-VCF-3514-S	VERTICALE CERVICAL FA Scr 3.5 x 14 mm solid	
S-VCF-3516-S	VERTICALE CERVICAL FA Scr 3.5 x 16 mm solid	
S-VCF-3518-S	VERTICALE CERVICAL FA Scr 3.5 x 18 mm solid	
S-VCF-3520-S	VERTICALE CERVICAL FA Scr 3.5 x 20 mm solid	
S-VCF-3522-S	VERTICALE CERVICAL FA Scr 3.5 x 22 mm solid	
S-VCF-3524-S	VERTICALE CERVICAL FA Scr 3.5 x 24 mm solid	
S-VCF-3526-S	VERTICALE CERVICAL FA Scr 3.5 x 26 mm solid	
S-VCF-3528-S	VERTICALE CERVICAL FA Scr 3.5 x 28 mm solid	
S-VCF-3530-S	VERTICALE CERVICAL FA Scr 3.5 x 30 mm solid	



System:
VERTICALE CERVICAL

Implant type:
Far angle screw

Typing:
smooth, cannulated

Material:
Ti6Al4V ELI

All articles are sterile
packed and include a set
screw.



Article number	Description	Figure
S-VCF-4020-KG	VERTICALE CERVICAL FA smooth 4 x 20 mm can	
S-VCF-4022-KG	VERTICALE CERVICAL FA smooth 4 x 22 mm can	
S-VCF-4024-KG	VERTICALE CERVICAL FA smooth 4 x 24 mm can	
S-VCF-4026-KG	VERTICALE CERVICAL FA smooth 4 x 26 mm can	
S-VCF-4028-KG	VERTICALE CERVICAL FA smooth 4 x 28 mm can	
S-VCF-4030-KG	VERTICALE CERVICAL FA smooth 4 x 30 mm can	
S-VCF-4032-KG	VERTICALE CERVICAL FA smooth 4 x 32 mm can	
S-VCF-4034-KG	VERTICALE CERVICAL FA smooth 4 x 34 mm can	
S-VCF-4036-KG	VERTICALE CERVICAL FA smooth 4 x 36 mm can	

VERTICALE® CERVICAL Implants

Article number	Description	Figure
S-VCF-4014-K	VERTICALE CERVICAL FA Scr 4.0 x 14 mm can	
S-VCF-4016-K	VERTICALE CERVICAL FA Scr 4.0 x 16 mm can	
S-VCF-4018-K	VERTICALE CERVICAL FA Scr 4.0 x 18 mm can	
S-VCF-4020-K	VERTICALE CERVICAL FA Scr 4.0 x 20 mm can	
S-VCF-4022-K	VERTICALE CERVICAL FA Scr 4.0 x 22 mm can	
S-VCF-4024-K	VERTICALE CERVICAL FA Scr 4.0 x 24 mm can	
S-VCF-4026-K	VERTICALE CERVICAL FA Scr 4.0 x 26 mm can	
S-VCF-4028-K	VERTICALE CERVICAL FA Scr 4.0 x 28 mm can	
S-VCF-4030-K	VERTICALE CERVICAL FA Scr 4.0 x 30 mm can	
S-VCF-4032-K	VERTICALE CERVICAL FA Scr 4.0 x 32 mm can	
S-VCF-4034-K	VERTICALE CERVICAL FA Scr 4.0 x 34 mm can	
S-VCF-4036-K	VERTICALE CERVICAL FA Scr 4.0 x 36 mm can	
S-VCF-4038-K	VERTICALE CERVICAL FA Scr 4.0 x 38 mm can	
S-VCF-4040-K	VERTICALE CERVICAL FA Scr 4.0 x 40 mm can	
S-VCF-4042-K	VERTICALE CERVICAL FA Scr 4.0 x 42 mm can	
S-VCF-4044-K	VERTICALE CERVICAL FA Scr 4.0 x 44 mm can	
S-VCF-4046-K	VERTICALE CERVICAL FA Scr 4.0 x 46 mm can	
S-VCF-4048-K	VERTICALE CERVICAL FA Scr 4.0 x 48 mm can	
S-VCF-4050-K	VERTICALE CERVICAL FA Scr 4.0 x 50 mm can	
S-VCF-4052-K	VERTICALE CERVICAL FA Scr 4.0 x 52 mm can	
S-VCF-4054-K	VERTICALE CERVICAL FA Scr 4.0 x 54 mm can	
S-VCF-4056-K	VERTICALE CERVICAL FA Scr 4.0 x 56 mm can	
S-VCF-4520-K	VERTICALE CERVICAL FA Scr 4.5 x 20 mm can	
S-VCF-4525-K	VERTICALE CERVICAL FA Scr 4.5 x 25 mm can	
S-VCF-4530-K	VERTICALE CERVICAL FA Scr 4.5 x 30 mm can	
S-VCF-4535-K	VERTICALE CERVICAL FA Scr 4.5 x 35 mm can	
S-VCF-4540-K	VERTICALE CERVICAL FA Scr 4.5 x 40 mm can	
S-VCF-4545-K	VERTICALE CERVICAL FA Scr 4.5 x 45 mm can	
S-VCF-4550-K	VERTICALE CERVICAL FA Scr 4.5 x 50 mm can	
S-VCF-4555-K	VERTICALE CERVICAL FA Scr 4.5 x 55 mm can	


System:
VERTICALE CERVICAL

Implant type:
Far angle screw

Typing:
cannulated

Material:
Ti6Al4V ELI

All articles are sterile packed and include a set screw.



VERTICALE® CERVICAL Implants








System:
VERTICALE CERVICAL

Implant type:
Set Screw / connector

Material:
Ti6Al4V ELI

All articles are sterile
packed.

* Package includes 2 pcs.

Article number	Description	Figure
S-VMS-1020	VERTICALE CERVICAL Set Screw X20*	
S-VCC-2833	VERTICALE CERVICAL Cross Connector 28-33 mm	
S-VCC-3241	VERTICALE CERVICAL Cross Connector 32-41 mm	
S-VCC-4057	VERTICALE CERVICAL Cross Connector 40-57 mm	
S-VCI-4035-I55	VERTICALE CERVICAL RC Inl 5.5 to 3.5/4 mm	
S-VCI-4035-T55	VERTICALE CERVICAL RC Tri 5.5 to 3.5/4 mm	
S-VCL-TUL	VERTICALE CERVICAL Lateral TUL Conn	

System:
VERTICALE CERVICAL

Implant type:
Rod

Typing:
Prebent, straight,
transition

Material:
Ti6Al4V ELI

All articles are sterile
packed.

All packages include 2 pcs.

Article number	Description	Figure
S-VOR-3520-P55	VERTICALE CERVICAL Rod Pre 3.5/200 mm 55°	
S-VOR-4020-P55	VERTICALE CERVICAL Rod Pre 4.0/200 mm 55°	
S-VCR-3506-S	VERTICALE CERVICAL Rod 3.5/60 mm	
S-VCR-3512-S	VERTICALE CERVICAL Rod 3.5/120 mm	
S-VCR-3524-S	VERTICALE CERVICAL Rod 3.5/240 mm	
S-VCR-4006-S	VERTICALE CERVICAL Rod 4.0/60 mm	
S-VCR-4012-S	VERTICALE CERVICAL Rod 4.0/120 mm	
S-VCR-4024-S	VERTICALE CERVICAL Rod 4.0/240 mm	
S-VCR-5535-T	VERTICALE CERVICAL Rod Trans. 5.5/3.5 mm	
S-VCR-5540-T	VERTICALE CERVICAL Rod Trans. 5.5/4.0 mm	

 D60035

VERTICAL[®] CERVICAL Instruments

Article number	Description	Figure
VI-0010	VERTICAL CERVICAL Awl with stop	
VI-0020	VERTICAL CERVICAL Probe, straight	
VI-0025	VERTICAL CERVICAL Pedicle Feeler	
VI-0030	VERTICAL CERVICAL Depth Gauge	
VI-0035	VERTICAL CERVICAL Tap 3.5 mm std	
VI-0040	VERTICAL CERVICAL Tap 4.0 mm std	
VI-0041	VERTICAL CERVICAL Tap 4.0 mm st+can	
VI-0130	VERTICAL CERVICAL Pedicle SD	
VI-0135	VERTICAL CERVICAL Tap 3.5 mm XL	
VI-0140	VERTICAL CERVICAL Tap 4.0 mm XL	
VI-0141	VERTICAL CERVICAL Tap 4.0 mm XL, can	
VI-0235	VERTICAL CERVICAL Drill 2.4 mm standard	
VI-0240	VERTICAL CERVICAL Drill 2.9 mm standard	
VI-0335	VERTICAL CERVICAL Drill 2.4 mm XL	
VI-0340	VERTICAL CERVICAL Drill 2.9 mm XL	

VERTICAL[®] CERVICAL Instruments

Article number	Description	Figure
VI-0201	VERTICAL CERVICAL Straight Handle, can	
VI-0211	VERTICAL CERVICAL Straight Handle Ratchet can	
VI-0230	VERTICAL CERVICAL Drill Guide Std	
VI-0256	VERTICAL CERVICAL Drill Guide XL	
VI-0258	VERTICAL CERVICAL Guide Wire 1.3 x 470 mm	
VI-0260	VERTICAL CERVICAL Rod Cutter	
VI-0261	VERTICAL CERVICAL Rod Cutter, in-situ	
VI-0292	VERTICAL CERVICAL Drill 2.4 mm NAV DG	
VI-0293	VERTICAL CERVICAL Drill 2.9 mm NAV DG	
VI-0294	VERTICAL CERVICAL Tap 3.5 mm NAV DG	
VI-0295	VERTICAL CERVICAL Tap 4.0 mm NAV DG	
VI-0296	VERTICAL CERVICAL Tap 4.0 mm NAV DG c	
VI-0301	VERTICAL CERVICAL T-Handle, cannulated	
VI-0311	VERTICAL CERVICAL T-Handle, Ratchet	
VI-0330	VERTICAL CERVICAL Rod Holder	
VI-0350	VERTICAL CERVICAL Rod and Tulip Adjuster	

VERTICALE® CERVICAL Instruments

Article number	Description	Figure
VI-0360	VERTICALE CERVICAL Reduction Instrument	
VI-0421	VERTICALE CERVICAL Setscrew Starter dbl	
VI-0440	VERTICALE CERVICAL Torque Limiter 3Nm	
VI-0446	VERTICALE CERVICAL Pedicle Screwdriver Ballhead	
VI-0450	VERTICALE CERVICAL Counter Torque	
VI-0535	VERTICALE CERVICAL Phantom Rod 3.5 x 200 mm	
VI-0610	VERTICALE CERVICAL Rod Bender	
VI-0611	VERTICALE CERVICAL Rod Bending Tube	
VI-0620	VERTICALE CERVICAL Distraction Pliers	
VI-0630	VERTICALE CERVICAL Compression Pliers	
VI-0700	VERTICALE CERVICAL Occiput Plate Holder	
VI-0710	VERTICALE CERVICAL OC Plate Bender right	
VI-0711	VERTICALE CERVICAL OC Plate Bender left	
VI-0720	VERTICALE CERVICAL OC Drill Guide	
VI-0730	VERTICALE CERVICAL OC Drill, straight	
VI-0731	VERTICALE CERVICAL OC Drill, flexible	

VERTICAL® CERVICAL Instruments

Article number	Description	Figure
VI-0740	VERTICAL CERVICAL OC Taps, straight	
VI-0741	VERTICAL CERVICAL OC Taps, flexible	
VI-0750	VERTICAL CERVICAL OC Screwdriver	
VI-0751	VERTICAL CERVICAL OC Screwdriver, flex.	
VI-0752	VERTICAL CERVICAL OC Screwdriver short	
VI-0760	VERTICAL CERVICAL OC Torque Limiter 3 Nm	
VI-0761	VERTICAL CERVICAL OC Counter Torque	
VI-0810	VERTICAL CERVICAL CC Torque Limiter 2Nm	
VI-0820	VERTICAL CERVICAL CC Counter Torque	
VI-0830	VERTICAL CERVICAL Cross Connector Sizer	
VI-0910	VERTICAL CERVICAL Nav Awl	
VI-0920	VERTICAL CERVICAL Nav Probe	
VI-0930	VERTICAL CERVICAL Nav Pedicle Screwdriver	
VI-0940	VERTICAL CERVICAL Nav Drill Guide	

Notes

Notes



 **Silony Medical GmbH**
Leinfelder Strasse 60
70771 Leinfelden-Echterdingen
Germany
Tel +49 711 78 25 25 0
Fax +49 711 78 25 25 11
www.silonyspine.com

 elabeling.silony-medical.com/contact

This information is only intended for the USA and its territories.

D60036.d.EN 10.02.2025

