

Instruction for Use (IFU050) - ALTOS® PCT Spinal Implants

ALTOS PCT

Important Information – please read before use!

These Instructions for Use apply to the ALTOS PCT device. Other sets of manufacturer’s instructions are available containing information required for application of the systems:

- Reprocessing Instructions for Instruments (LBL379)
- Instrumentation Guide (STG018)

The system-specific product information can also be viewed at the following links:

- <https://elabeling.silony-medical.com>
- www.silonyspine.com

Table of Contents

1	Purpose	1
2	Device Description.....	1
3	Indications.....	2
4	Contraindications	2
5	Warnings and Precautions	2
6	Potential adverse events with ALTOS PCT	3
5	Implant Selection	4
5.1	Pre-Operative	4
5.2	Intra-Operative	4
5.3	Post-Operative	5
6	Packaging.....	6
7	Recommendations for Steam Sterilization	6
8	Cleaning and Decontamination.....	6
9	Magnetic Resonance Imaging (MRI) Safety Information	7
10	Glossary of Symbols	7

1 Purpose

ALTOS® Posterior Cervical Thoracic (PCT) System is intended to help provide immobilization and stabilization of spinal segments as an adjunct to fusion of the cervical and/or upper thoracic spine (C1-T3).

2 Device Description

ALTOS PCT is a posterior cervical-thoracic system which consists of a variety of shapes and sizes of rods, hooks, polyaxial screws, and connecting components, which can be rigidly locked to the rod in a variety of configurations. Each construct can be tailor-made for the individual case.

The ALTOS PCT system is fabricated from medical grade titanium alloy (Ti6Al4V) according to ASTM F136 or ISO 5832-3.

3 Indications

ALTOS Posterior Cervical Thoracic System is intended to provide immobilization and stabilization of spinal segments as an adjunct to fusion for the following acute and chronic instabilities of the cervical spine (C1 to C7) and the thoracic spine from T1-T3: traumatic spinal fractures and/or traumatic dislocations; instability or deformity; failed previous fusions (e.g., pseudarthrosis); tumors involving the cervical spine; and degenerative disease, including intractable radiculopathy and/or myelopathy, neck and/or arm pain of discogenic origin as confirmed by radiographic studies, and degenerative disease of the facets with instability. The ALTOS system is also intended to restore the integrity of the spinal column even in the absence of fusion for a limited time period in patients with advanced stage tumors involving the cervical spine in whom life expectancy is of insufficient duration to permit achievement of fusion.

4 Contraindications

Contraindications include, but are not limited to:

- Active infectious process or significant risk of infection (immunocompromised).
- Signs of local inflammation.
- Fever or leukocytosis.
- Morbid obesity.
- Pregnancy.
- Mental illness.
- Grossly distorted anatomy caused by congenital abnormalities.
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count.
- Suspected or documented metal allergy or intolerance.
- Any case not needing a bone graft and fusion.
- Any case where the implant components selected for use would be too large or too small to achieve a successful result.
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality.
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.
- Any patient unwilling to follow postoperative instructions.
- Any case not described in the indications.

NOTA BENE: Although not absolute contraindications, conditions to be considered as potential factors for not using this device include:

- Severe bone resorption.
- Osteomalacia.
- Severe osteoporosis.

5 Warnings and Precautions

- The implantation of spinal fixation systems should be performed only by experienced spinal surgeons with specific training in the use of these spinal systems because this is a technically demanding procedure presenting a risk of serious injury to the patient. Preoperative planning and patient anatomy should be considered when selecting implant diameter and length.
- Patients with previous spinal surgery at the levels to be treated may not experience the same clinical outcomes as those without a previous surgery.
- Selection of an appropriately sized device for the patient is important and increases the likelihood of a satisfactory outcome.

- The implantation of screw based spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this type of device.
- Instrumentation provided with the implants must be used in accordance with the approved surgical technique.
- Do not use excessive force when introducing and positioning the implant to avoid damaging the implant.
- Re-usable surgical instruments must be re-sterilized prior to next use.
- Do not reuse the device even if the device shows no external signs of damage. Internal stresses from previous use may cause early failure.
- Should not be used with components of any other system or manufacturer.
- Do not use titanium with stainless steel in the same construct. Premature device failure and / or infection in the patient may occur.
- Physicians should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may have an impact on the performance of this system.
- This device should be used in conjunction with a bone graft.
- The ALTOS system has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safe use of the ALTOS system in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

6 Potential adverse events with ALTOS PCT

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of potential adverse events includes, but is not limited to:

- Early or late loosening of any or all of the components.
- Disassembly, bending, and/or breakage of any or all of the components.
- Foreign body (allergic) reaction to implants, debris, corrosion products (from crevice, fretting, and/or general corrosion), including metallosis, staining, tumor formation, and/or autoimmune disease.
- Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, fibrosis, necrosis, and/or pain. Bursitis. Tissue or nerve damage caused by improper positioning and placement of implants or instruments.
- Post-operative change in spinal curvature, loss of correction, height, and/or reduction.
- Infection.
- Dural tears, pseudomeningocele, fistula, persistent CSF leakage, meningitis.
- Loss of neurological function (e.g., sensory and/or motor), including paralysis (complete or incomplete), dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits.
- Neuropathy, neurological deficits (transient or permanent), paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss.
- Urinary retention or loss of bladder control or other types of urological system compromise.
- Scar formation possibly causing neurological compromise or compression around nerves and/or pain.
- Fracture, microfracture, resorption, damage, or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above, and/or below the level of surgery.
- Retropulsed graft.
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery.
- Non-union (or pseudarthrosis). Delayed union. Mal-union.
- Loss of or increase in spinal mobility or function.
- Inability to perform the activities of daily living.
- Bone loss or decrease in bone density, possibly caused by stresses shielding.

- Graft donor site complications including pain, fracture, or wound healing problems.
- Ileus, gastritis, bowel obstruction or loss of bowel control or other types of gastrointestinal system compromise.
- Hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise.
- Reproductive system compromise, including sterility, loss of consortium, and sexual dysfunction.
- Development of respiratory problems, e.g., pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.
- Change in mental status.
- Death.

NOTE: Additional surgery may be necessary to correct some of these potential adverse events.

5 Implant Selection

The selection of the proper size, shape and design of the implant for each patient is crucial to the success of the procedure. Metallic surgical implants are subject to repeated stresses in use, and their strength is limited by the need to adapt the design to the size and shape of human bones. Unless great care is taken in patient selection, proper placement of the implant, and postoperative management to minimize stresses on the implant, such stresses may cause metal fatigue and consequent breakage, bending or loosening of the device before the healing process is complete, which may result in further injury or the need to remove the device prematurely.

5.1 Pre-Operative

- Only patients that meet the criteria described in the indications should be selected.
- Patient conditions and/or pre-dispositions such as those addressed in the aforementioned contraindications should be avoided.
- Care should be used in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Implants and instruments should be protected during storage, especially from corrosive environments.
- An adequate inventory of implants should be available at the time of surgery, normally a quantity in excess of what is expected to be used.
- Since mechanical parts are involved, the surgeon should be familiar with the various components before using the equipment and should personally assemble the devices to verify that all parts and necessary instruments are present before the surgery begins. ALTOS Posterior Cervical Thoracic System components (described in the Device Description section) are not to be combined with the components from another manufacturer. Different metal types should never be used together.
- All components and instruments should be cleaned and sterilized before use. (See directions below). Additional sterile components should be available in case of an unexpected need.

5.2 Intra-Operative

- The surgeon follows the surgical technique and instructions for use of the device.
- All components should be inspected and determined to be free of damage.
- Extreme caution should be used around the spinal cord and nerve roots. Damage to the nerves will cause loss of neurological functions.
- Breakage, slippage, or misuse of instruments or implant components may cause injury to the patient or operative personnel.
- The rods should not be repeatedly or excessively bent. The rods should not be reverse bent in the same location. Use great care to ensure that the implant surfaces are not scratched or notched, since such actions may reduce the functional strength of the construct. If the rods are cut to length, they should be cut in such a way as to create a

flat, non-sharp surface perpendicular to the midline of the rod. Cut the rods outside the operative field. Whenever possible, use pre-cut rods of the length needed.

- Whenever possible or necessary, an imaging system should be utilized to facilitate surgery.
- To insert a screw properly, drill a pilot hole corresponding to selected screw size and prepare screw site with a sharp tap.
- **Caution:** Do not overlap or use a screw that is either too long or too large. Overtapping or using an incorrectly sized screw may cause nerve damage, hemorrhage, or the other possible adverse events listed elsewhere in this package insert.
- Before closing the soft tissues, all of the screws or set screws should be tightened firmly. Recheck the tightness of all screws or set screws after finishing to make sure that none loosened during the tightening of the other screws or set screws. Failure to do so may cause loosening of the other components.

5.3 Post-Operative

The physician's postoperative directions and warnings to the patient, and the corresponding patient compliance, are extremely important.

- The choice to administer post-operative antibiotics is at the discretion of the surgeon.
- Post-operative mobilization and rehabilitation is at the discretion of the surgeon dependent on clinical and radiological progress.
- The need for external orthotic support is not mandatory with the final choice based on surgeon preference, patient condition and intra-operative findings that might influence implant security.
- The patient is to be instructed to reduce undue stress on the implant as a precaution to avoid clinical problems that could result in fixation failure.
- The patient is to be instructed to follow the post-operative regime.
- Optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The MR image quality of artifacts may be compromised if the area of interest is in the exact same area or relatively close to the position of the ALTOS implant.
- ALTOS Posterior Cervical Thoracic System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safe use of the ALTOS system in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
- Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening and/or breakage of the device(s) are complications which may occur as a result of excessive or early weight-bearing or muscular activity. The risk of bending, loosening, or breakage of a temporary internal fixation device during postoperative rehabilitation may be increased if the patient is active, or if the patient is debilitated or demented. The patient should be warned to avoid falls or sudden jolts in spinal position.
- To allow the maximum chances for a successful surgical result, the patient or devices should not be exposed to mechanical vibrations or shock that may loosen the device construct. The patient should be warned of this possibility and instructed to limit and restrict physical activities, especially lifting and twisting motions and any type of sport participation. The patient should be advised not to smoke tobacco or utilize nicotine products, or to consume alcohol or non-steroidal or anti-inflammatory medications such as aspirin during the bone graft healing process.
- The patient should be advised of their inability to bend or rotate at the point of spinal fusion and taught to compensate for this permanent physical restriction in body motion.
- Failure to immobilize a delayed or non-union of bone will result in excessive and repeated stresses on the implant. By the mechanism of fatigue, these stresses can cause the eventual bending, loosening, or breakage of the device(s). It is important that immobilization of the spinal surgical site be maintained until firm bony union is established and confirmed by roentgenographic examination. If a state of non-union persists or if the components loosen, bend, and/or break, the device(s) should be

revised and/or removed immediately before serious injury occurs. The patient must be adequately warned of these hazards and closely supervised to ensure cooperation until bony union is confirmed.

- As a precaution, before patients with implants receive any subsequent surgery (such as dental procedures), prophylactic antibiotics may be considered, especially for high-risk patients.
- The ALTOS Posterior Cervical Thoracic System implants are temporary internal fixation devices. Internal fixation devices are designed to stabilize the operative site during the normal healing process. After the spine is fused, these devices serve no functional purpose and should be removed. While the final decision on implant removal is, of course, up to the surgeon and patient, in most patients, removal is indicated because the implants are not intended to transfer or support forces developed during normal activities. If the device is not removed following completion of its intended use, one or more of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position, possibly resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening and breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; (7) Bone loss due to stress shielding; and (8) Potential unknown and/ or unexpected long term effects such as carcinogenesis. Implant removal should be followed by adequate postoperative management to avoid fracture, re-fracture, or other complications.
- Any retrieved devices should be treated in such a manner that reuse in another surgical procedure is not possible. As with all orthopedic implants, the ALTOS Posterior Cervical Thoracic System components should never be reused under any circumstances.

6 Packaging

Packages for each of the components should be intact upon receipt. If a loaner or consignment system is used, all sets should be carefully checked for completeness and all components including instruments should be carefully checked to ensure that there is no damage prior to use. Damaged packages or products should not be used and should be returned to Silony Spine.

7 Recommendations for Steam Sterilization

All implants of the ALTOS Posterior Cervical Thoracic System are provided non-sterile and must be cleaned and sterilized prior to uses. Remove all packaging materials prior to sterilization. Only FDA-cleared sterilization packaging materials should be used when sterilizing the ALTOS Posterior Cervical Thoracic System implants and instruments. Only sterile products should be placed in the operative field. Unless specified elsewhere, these products are recommended to be steam sterilized by the hospital using the process parameters below:

METHOD	CYCLE	TEMPERATURE	EXPOSURE	DRY TIMES
Steam	Pre-Vacuum	270°F (132°C)	4 Minutes	30 Minutes

NOTE: Because of the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g. temperatures, times) used for their equipment.

8 Cleaning and Decontamination

Unless just removed from an unopened sterile Silony Spine package, all instruments and implants must be unpackaged, disassembled (if applicable) and cleaned before sterilization and introduction into a sterile surgical field or, if applicable, return of the product to Silony Spine. Remove all packaging materials prior to disassembly (if applicable) and cleaning. Cleaning instructions and

associated disassembly instructions (if applicable) can be found in LBL379, "Care, Maintenance and Sterilization of Non-Sterile Implants and Instruments."











NOTE: Certain cleaning solutions such as those containing formalin, glutaraldehyde, bleach and/or other alkaline cleaners may damage some devices, particularly instruments; these solutions should not be used unless otherwise instructed by Silony Spine.

9 Magnetic Resonance Imaging (MRI) Safety Information


ALTOS Posterior Cervical Thoracic System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safe use of the ALTOS system in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

10 Glossary of Symbols

Symbols according to ISO 15223-1 (Recognition number: 5-134) and 21 CFR 801.109:

Symbol	Explanatory Text
	Manufacturer
	Catalog number
	Batch code
	Consult instructions for use
	Caution
	Do not reuse
	Sterilized using irradiation
	Federal law in the USA restricts this device to sale by or on the order of a physician
	Contact
	Caution MR* uncertain

*The term MR is synonymous with MRI and means magnetic resonance imaging

 Silony Spine Corp.
 8200 NW 27th St
 STE#104
 Doral, FL 33122
 Tel: +1 305 916 0016
www.silonyspine.com