

Resection of the facet joint

SI-0019 SPECIALS Chisel short 10 mm width



The transforminal approach to the disc space is usually achieved by means of unilateral resection of the facet joint on the approach side (Fig. 1). The 10 mm SPECIALS Chisel short can be used for this, as well as standard instruments such as Luer forceps and rongeurs.




Fig. 1 Resection of a facet joint with the short chisel

D30003

⚠ D30001

D30111.a.EN 21.03.2018

Silony Medical Europe GmbH

Bahnhofstraße 1 28195 Bremen

Tel +49 421 246 956 0 Fax +49 421 246 956 55

Silony Medical GmbH

Leinfelder Straße 60 70771 Leinfelden-Echterdingen Tel +49 711 78 25 25 0

Fax +49 711 78 25 25 0

